



MEDICAL INC.

2692 HUNTLEY ■ MEMPHIS, TN 38132 ■ (901) 345-0606 ■ FAX (901) 345-0608  
E-MAIL [sales@kenad.com](mailto:sales@kenad.com) [www.kenad.com](http://www.kenad.com)

January 2000

Dear Valued Customer,

It is with enthusiasm that Kenad presents our most comprehensive catalog to date. You will find the new user-friendly format most helpful in locating the products you need.

One call gives you access to the widest selection of orthopedic surgical support and patient protective devices. Our knowledgeable and experienced staff is always available to answer your questions or help you choose the right products. At Kenad Medical, whether a large or small company, you can benefit from our commitment to priority service and to the individual attention we give each customer. You can always depend on Kenad Medical; we are here to serve you.

As always, your suggestions for modifications, improvements, and new ideas are welcome. We realize you have a choice in which company you choose. We appreciate Kenad Medical as your choice.

Cordially Yours,

Kenneth Elmore  
President

THE SOFTGOODS SPECIALIST

# GENERAL INFORMATION

## Obtaining Special Products ►

Products listed in the catalog are stock items. Special products may be manufactured upon request. Please call for ordering and pricing information. Special ordered merchandise may not be returned. Please allow 10 working days for delivery of special order merchandise.

## Minimum Order ►

There is a minimum of \$25.00 per order. Orders less than the minimum order amount will be charged a \$5.00 processing fee.

## Terms ►

The following terms and conditions are intended to outline some good business practice. They are normal practices and purchasing terms; your adherence to these terms and practices will allow us to continue to offer excellent service at prices well below the prices of our major competitors: Net 30 days. There will be a one percent discount for payment within 15 days of invoice date. Accounts more than 45 days past due will be placed on credit hold; accounts more than 60 days past due will be placed on C.O.D. status. New accounts will ship C.O.D. pending credit approval. A \$2.00 handling charge will be added to each order.

## Shipping ►

All prices are shown F.O.B. Memphis, Tennessee. Orders will be shipped prepaid with the shipping charge added to the invoice. There is a \$10.00 charge for drop shipments.

## Examination of Merchandise ►

Upon receipt of the order, the customer should examine it carefully to be certain that sizes, quantities, etc., agree with the packing list. CLAIMS FOR SHORTAGES MUST BE MADE WITHIN FIVE DAYS AFTER RECEIPT OF ORDER.

## Limited Warranty ►

All KENAD products are warranted for 90 days from date of invoice as to product quality and workman ship. Any defective merchandise will, at the option of KENAD, be repaired or replaced.

Any implied warranties, including warranties of merchantability or fitness for a particular purpose, which arise out of the sale of a KENAD product are limited to a period of 90 days from date of invoice. These written warranties are given in lieu of any implied warranties, including warranties of merchantability or fitness for a particular purpose.

## Prices ►

Prices are subject to change without notice.

## Specifications ►

Specifications may be changed to comply with more current standards.

## Return Goods Policy ►

Merchandise may be returned within 60 days following receipt of order. Nothing may be returned without prior authorizations. Return authorization number must appear on all packages being returned, Copies of the original invoice must accompany the return. Packages returned without a return authorization number will be refused. Restocking charges for returns other than those which are the result of our error will be charged as follows:

One to 60 days after customer receipt of order – 25%

**\*\*RETURNS WHICH ARE MORE THAN 60 DAYS OLD WILL NOT BE ACCEPTED.**

# GENERAL FITTING GUIDE

This chart provides general sizing guidelines only. The patients' individual physical characteristics may require a different size product.

## Knee Braces & Supports

*Measure around center of knee with leg extended.*

XSmall	12" - 13"
Small	13" - 14"
Medium	14" - 15"
Large	15" - 16"
XLarge	16" - 18"
XXLarge	18" - 20"

## Wrist & Forearm

*Measure around wrist.*

Small	5" - 6"
Medium	6" - 7"
Large	7" - 8"

## Ankle Support

*By shoe size.*

XSmall	5 - 6
Small	6 - 7
Medium	8 - 9
Large	10 - 12
XLarge	12 - 14

*Size may vary according to ankle circumference.*

## Lumbar Support

*Measure around waist.*

XSmall	22" - 26"
Small	26" - 30"
Medium	30" - 34"
Large	34" - 38"
XLarge	38" - 42"
XXLarge	42" - 46"

## Elbow Support

*Measure the largest circumference of the forearm.*

XSmall	8" - 9"
Small	9" - 10"
Medium	10" - 11"
Large	11" - 12"
XLarge	12" - 14"

## Short Leg Walker Brace

*Equivalent Shoe Size.*

	Men	Women
Small	6 - 8	6 - 8½
Medium / Large	8½ - 10½	9 - 11
XLarge	11 - 13	11½ - 13

*Fits left or right.*



The charts provide general sizing guidelines only – the patient's individual physical characteristics may require a different size product. The proper medical authority must make the final decision on what size product to use.

**Patient Protective Vest Adult Sizing Chart**

Binding Color	Size	Patient Weight Lb.	Weight Kg.	Chest In.	Size Cm.
White	XSmall	60 - 90	27 - 41	28 - 34	71 - 86
Red	Small	80 - 120	36 - 54	32 - 39	81 - 99
Green	Medium	120 - 155	50 - 70	37 - 45	94 - 114
Yellow	Large	145 - 190	66 - 86	43 - 49	109 - 124
Blue	XLarge	180 - 230	62 - 104	47 - 52	119 - 132
Black	XXLarge	100 - 125	100 - 125	50 - 56	127 - 142

**Toddler Sizing Chart**

Color	Size	Age
Red	Small	2 - 4 Years
Green	Medium	4 - 6 Years
Blue	Large	6 - 8 Years

# TABLE OF CONTENTS

	 <u>Page</u>
Knee Immobilizer .....	1-3
Knee Braces & Supports .....	4-10
Ankle Supports .....	11-16
Cast Room Products .....	17-25
Splints & Supports .....	26-28
Wrist Splints & Supports .....	29-37
Kenad Rib Belts .....	38
Kenad Abdominal Binders .....	39
Kenad Lumbosacral Supports .....	40-45
Clavicular Supports .....	46-50
Cervical Collars .....	51-53
Kenad Shoulder Immobilizers, Slings .....	54-56
Kenad Arm Slings .....	57-62
Kenad Traction Products .....	63-66
Patient Protective Products .....	67-72
Patient Assist Products .....	73-79
Splints .....	80-85
Rehabilitation .....	86



# KNEE IMMOBILIZERS

Kenad SG Medical, Inc.

Knee immobilizers may be used postoperatively, after arthroscopic procedures, after cast removal, or in place of a cast, and in conservative treatment of sprains, strains, ligamentous injuries, disease, or deteriorating process.

## Universal Knee Immobilizer



Size	Catalog No.
12"	14-1012
16"	14-1016
19"	14-1019
20"	14-1020
22"	14-1022
24"	14-1024
26"	14-1026
28"	14-1028

### Construction

$\frac{3}{8}$ " foam with hook compatible soft brushed pile nylon. Spring steel stays with loop lock closure.

### Features

May be trimmed for custom fit. Five sturdy preformed spring steel stays for added stability.

## Universal Knee Immobilizer w/Elastic Closure



Size	Catalog No.
12"	14-1112
16"	14-1116
19"	14-1119
23"	14-1123

### Construction

Soft blue foam, spring steel stays with elastic closures.

### Features

Elastic allows for added compression. Spring steel stays for extra rigidity.

## Economy Universal Knee Immobilizer



Size	Catalog No.
12"	14-1212
16"	14-1216
19"	14-1219
23"	14-1223

### Construction

$\frac{3}{8}$ " foam with hook compatible soft brushed pile nylon. Preformed aluminum stays with loop lock closure.

### Features

Movable stays that can be adjusted to fit patient. foam may be trimmed for sized custom fit. Seven aluminum stays for added stability.

## Deluxe 3-Panel Universal Knee Immobilizer



Size	Catalog No.
12"	14-2012
16"	14-2016
19"	14-2019
22"	14-2022
24"	14-2024

### Construction

Center panel made of  $\frac{3}{8}$ " foam with hook compatible soft brushed pile nylon. Side panel made of canvas foam laminate. Preformed spring steel stays. Hook and loop closure.

### Features

Adjustable side panels may be adjusted to allow full range of sizes. Foam may be trimmed for sized custom fit. Five sturdy spring steel stays for added stability.

# KNEE IMMOBILIZERS

Kenod SG Medical, Inc.

## Universal Cylinder Splint



Size	Catalog No.
16"	14-3016
19"	14-3019
22"	14-3022
24"	14-3024

### Construction

Blue velfoam with movable hook and loop attached cotton webbing straps with adjustable loop-lock closures. Nine galvanized steel stays. Pre-bent proximal, lateral, and anterior stays.

### Features

Commonly applied from proximal thigh to just above the ankle for knee injuries. May be used in place of cylinder cast. Nine sturdy spring steel stays for super rigid immobilization.

## Deluxe Knee Immobilizer



	12" Length	16" Length	19" Length	22" Length	24" Length
Size	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.
XSmall	15-1003	15-1503	15-2003	15-4003	15-4503
Small	15-1004	15-1504	15-2004	15-4004	15-4504
Medium	15-1005	15-1505	15-2005	15-4005	15-4505
Large	15-1006	15-1506	15-2006	15-4006	15-4506
XLarge	15-1007	15-1507	15-2007	15-4007	15-4507

### Construction

$\frac{3}{8}$ " canvas foam laminate. Preformed spring steel stays. Hook and loop closure.

### Features

Sturdy spring steel stays assure immobilization when properly applied. Available in multiple lengths and sizes that permit specific fitting of individual patients.

## Felt-Lined Knee Immobilizer



	12" Length	16" Length	19" Length	23" Length
Size	Catalog No.	Catalog No.	Catalog No.	Catalog No.
XSmall	15-2103	15-2503	15-2703	15-2603
Small	15-2104	15-2504	15-2704	15-2604
Medium	15-2105	15-2505	15-2705	15-2605
Large	15-2106	15-2506	15-2706	15-2606
XLarge	15-2107	15-2507	15-2707	15-2607

### Construction

$\frac{3}{8}$ " canvas foam laminate. Lined with  $\frac{1}{4}$ " orthopedic felt. Preformed spring steel stays. Hook and loop closure.

### Features

Sturdy spring steel stays assure immobilization when properly applied. Felt provides added rigidity. Available in multiple lengths and sizes that permit specific fitting of individual patients.

# KNEE IMMOBILIZERS

Kenad SG Medical, Inc.

## Economy 3-Panel Universal Knee Immobilizer



Size	Catalog No.
12"	14-6000
16"	14-6001
20"	14-6002
24"	14-6003

### Construction

Center panel made of  $\frac{3}{8}$ " foam with hook compatible soft brushed pile nylon. side panel made of canvas foam laminate. Preformed aluminum stays.

### Features

Adjustable side panels may be adjusted to allow full range of sizes. Foam may be trimmed for sized custom fit. Five sturdy aluminum stays for added stability. Hook and loop closure.

## Perforated Vinyl Knee Immobilizer



Size	12" Length Catalog No.	16" Length Catalog No.	19" Length Catalog No.	23" Length Catalog No.
XSmall	15-2803	15-2903	15-3003	15-3103
Small	15-2804	15-2904	15-3004	15-3104
Medium	15-2805	15-2905	15-3005	15-3105
Large	15-2806	15-2906	15-3006	15-3106
XLarge	15-2807	15-2907	15-3007	15-3107

### Construction

$\frac{3}{8}$ " foam laminated to perforated vinyl. Spring steel stays with hook and loop closure.

### Features

Sturdy spring steel stays assure immobilization when properly applied. Available in multiple lengths and sizes that permit specific fitting of individual patients. Vinyl is easily cleaned with damp cloth.

## Economy Knee Immobilizer



Size	12" Length Catalog No.	16" Length Catalog No.	19" Length Catalog No.	23" Length Catalog No.
XSmall	15-5503	15-5603	15-5703	15-5803
Small	15-5504	15-5604	15-5704	15-5804
Medium	15-5505	15-5605	15-5705	15-5805
Large	15-5506	15-5606	15-5706	15-5806
XLarge	15-5507	15-5607	15-5707	15-5807

### Construction

$\frac{3}{8}$ " foam with hook compatible soft brushed pile nylon. Preformed aluminum stays. Hook and loop closure.

### Features

Sturdy spring steel stays assure immobilization when properly applied. Available in multiple lengths and sizes that permit specific fitting of individual patients.



# KNEE BRACES & SUPPORTS

Kenod SG Medical, Inc.

Used preoperatively or postoperatively to help support body weight and to limit joint motion while providing medial-lateral stability, and in cases of knee instability due to sprains, strains, ligamentous injuries, disease, and the like.

## Rom Post-Op Knee Brace



Size	Catalog No.
Long	16-0100L
Short	16-0100S

### Construction

Soft black foam. Hook and Loop closure. Segmented hinge arms.

### Features

No migration down the leg. Hinge can be cut with cast saw for custom fit. Controlled range of motion with easy-to-adjust pin stops.

## Osteo-Arthritis Brace



Size	Catalog No.
Left Small/Medium	16-0104
Left Large/XLarge	16-0106
Right Small/Medium	16-0204
Right Large/XLarge	16-0206

### Construction

4mm high tensile hinge arms with a standard four-degree valgus force applied.

### Features

Low profile bicentric hinge. Simple strapping maintains force during movement. Anatomical shell design for maximum positional stability.

## Neoprene Knee Support w/Tubular Buttress



Size	Catalog No.
XSmall	16-0603
Small	16-0604
Medium	16-0605
Large	16-0606
XLarge	16-0607
XXLarge	16-0608
XXXLarge	16-0609

### Construction

1/8" neoprene nylon two sides. Tubular neoprene donut.

### Features

Open patella and tubular neoprene donut buttress helps maintain patella femoral groove alignment. Ideal for minor strains, sprains, or after post cast removal and/or bracing. Provides warmth, support, and equal compression.

## Neoprene Patella Stabilizer w/Stays & Tubular Buttress



Size	Catalog No.
XSmall	16-0703
Small	16-0704
Medium	16-0705
Large	16-0706
XLarge	16-0707
XXLarge	16-0708
XXXLarge	16-0709

### Construction

1/8" neoprene nylon two sides. Tubular neoprene donut. Two medial and lateral spiral steel stays.

### Features

Open patella and tubular neoprene donut buttress helps maintain patella femoral groove alignment. Ideal for minor strains, sprains, or after post cast removal and bracing. Provides warmth, support, and equal compression.



# KNEE BRACES & SUPPORTS

Kenad SG Medical, Inc.

## 13" Challenger Knee Brace



Size	Catalog No.
XSmall	16-0803
Small	16-0804
Medium	16-0805
Large	16-0806
XLarge	16-0807
XXLarge	16-0808
XXXLarge	16-0809

### Construction

Constructed of  $\frac{3}{16}$ " neoprene nylon two sides. Medial and lateral high tech hinges with easy to adjust pin stops. Hook and loop strapping for added stability.

### Features

Open patella helps maintain patella femoral groove alignment. High tech hinges and hook and loop strapping provides greater stability. Controlled range of motion with easy to adjust pin stops. Provides warmth, support, and equal compression. Ideal for minor strains, sprains, after post-op surgery, cast removal, and bracing.

## 17" Challenger Knee Brace



Size	Catalog No.
XSmall	16-0903
Small	16-0904
Medium	16-0905
Large	16-1906
XLarge	16-0907
XXLarge	16-0908
XXXLarge	16-0909

### Construction

Constructed of  $\frac{1}{8}$ " neoprene nylon, two sides. Medial and lateral high tech hinges with easy to adjust pin stops. Hook and loop strapping for added stability.

### Features

Extra length and strapping provides a longer fulcrum for greater support. Open patella helps maintain patella femoral groove alignment. High tech hinges and hook and loop strapping provides greater stability. Controlled range of motion with easy to adjust pin stops. Provides warmth, support, and equal compression. Ideal for minor strains, sprains, after post-op surgery, cast removal, and bracing.

## 13" Challenger Knee Brace w/Anterior Closure



Size	Catalog No.
XSmall	16-1003
Small	16-1004
Medium	16-1005
Large	16-1006
XLarge	16-1007
XXLarge	16-1008
XXXLarge	16-1009

### Construction

Constructed of  $\frac{3}{16}$ " neoprene nylon, two sides. Medial and lateral high tech hinges with easy to adjust pin stops. Hook and loop strapping for added stability.

### Features

Anterior closure allows for easy on and easy off patients. Open patella helps maintain patella femoral groove alignment. High tech hinges and hook and loop strapping provides greater stability. Controlled range of motion with easy to adjust pin stops. Provides warmth, support, and equal compression. Ideal for minor strains, sprains, after post-op surgery, cast removal, and bracing.

# KNEE BRACES & SUPPORTS

Kenad SG Medical, Inc.

## Neoprene Knee Sleeve – Closed Patella



Size	Catalog No.
XSmall	16-1103
Small	16-1104
Medium	16-1105
Large	16-1106
XLarge	16-1107
XXLarge	16-1108
XXXLarge	16-1109

### Construction

1/8" neoprene nylon, two sides.

### Features

Provides warmth, support, and equal compression.

## Neoprene Knee Sleeve – Open Patella



Size	Catalog No.
XSmall	16-1203
Small	16-1204
Medium	16-1205
Large	16-1206
XLarge	16-1207
XXLarge	16-1208
XXXLarge	16-1209

### Construction

1/8" neoprene nylon, two sides.

### Features

Provides warmth, support, and equal compression.  
Patella cut out.

## Neoprene Patella Stabilizer w/Universal Buttress



Size	Catalog No.
XSmall	16-1303
Small	16-1304
Medium	16-1305
Large	16-1306
XLarge	16-1307
XXLarge	16-1308
XXXLarge	16-1309

### Construction

1/4" neoprene nylon two sides with horseshoe felt buttress. Two medial and lateral spiral steel stays Hook and loop strapping.

### Features

Provides warmth, support, and equal compression.  
Patella cut out. Horseshoe felt buttress can be positioned around the inside of the patella.

## Elastic Slip-on Knee Support



Size	Catalog No.
Small	16-1504
Medium	16-1505
Large	16-1506
XLarge	16-1507

### Construction

Heavy, two-ply elastic nylon/rubber blend for maximum conformability.

### Features

A slip-on knee support. Available in four sizes. Suntan color. 12" long.



# KNEE BRACES & SUPPORTS

Kenad SG Medical, Inc.

## Neoprene Hinged Knee Brace



Size	Catalog No.
XSmall	16-1603
Small	16-1604
Medium	16-1605
Large	16-1606
XLarge	16-1607
XXLarge	16-1608
XXXLarge	16-1609

### Construction

1/4" neoprene nylon, two sides. Horseshoe felt buttress and simple hinge.

### Features

Provides warmth, support, and equal compression. Universal felt buttress. Hinged for extra stability.

## Neoprene Thigh Sleeve



Size	Catalog No.
XSmall	16-1703
Small	16-1704
Medium	16-1705
Large	16-1706
XLarge	16-1707
XXLarge	16-1708
XXXLarge	16-1709

### Construction

1/8" neoprene nylon, two sides.

### Features

Provides warmth, support, and equal compression. Ideal for minor strains, sprains, after post-op surgery, case removal, and bracing.

## Elastic Cartilage Knee Brace



Size	Catalog No.
Small	16-2004
Medium	16-2005
Large	16-2006
XLarge	16-2007

### Construction

Heavy, two-ply elastic, nylon, and rubber blend. Medial and lateral stays are made of spiral spring steel to allow knee flexion. Removable cartilage pads give added support to patella.

### Features

Ideal knee support following injury or damage to cartilage, ligaments, and menisci.

## Neoprene Sleeve w/Lateral "J" Buttress



Size	Right Catalog No.	Left Catalog No.
XSmall	16-2103	16-2203
Small	16-2104	16-2204
Medium	16-2105	16-2205
Large	16-2106	16-2206
XLarge	16-2107	16-2207
XXLarge	16-2108	16-2208
XXXLarge	16-2109	16-2209

### Construction

1/8" neoprene nylon, two sides. Lateral tubular "J" buttress. Two medial and lateral spiral steel stays.

### Features

Provides warmth, support, and equal compression. Lateral tubular "J" buttress helps control femoral groove tracking of the patella.

# KNEE BRACES & SUPPORTS

Kenod SG Medical, Inc.

## Body Flex™ Knee Sleeve – Closed Patella



Size	Catalog No.
Small	16-2304
Medium	16-2305
Large	16-2306
XLarge	16-2307
XXLarge	16-2308

### Construction

2.1mm thick nylon lycra. Closed patella.

### Features

Less heat, breathable and thin as skin. Greater compression and non-slip suspension system.

## Body Flex™ Knee Sleeve Open Patella



Size	Catalog No.
Small	16-2404
Medium	16-2405
Large	16-2406
XLarge	16-2407
XXLarge	16-2408

### Construction

2.1mm thick nylon lycra. Open patella.

### Features

Less heat, breathable and thin as skin. Greater compression and non-slip suspension system.

## Body Flex™ Patella Brace w/Gel Buttfress



Size	Hinged	Stays
Small	16-2304H	16-2404S
Medium	16-2305H	16-2405S
Large	16-2306H	16-2406S
XLarge	16-2307H	16-2407S
XXLarge	16-2308H	16-2408S

### Construction

2.1mm thick nylon lycra with gel buttress. Available with hinges or spiral steel stays.

### Features

Less heat, breathable, and thin as skin. Greater compression and non-slip suspension system. Gel buttress helps control patella femoral groove. Tracking of the patella. Hinges or spiral steel stays for added stability.

## Hinged Elastic Knee Brace



Size	Catalog No.
Small	16-3004
Medium	16-3005
Large	16-3006
XLarge	16-3007

### Construction

Heavy, two-ply elastic and top-grade leather. Medial and lateral hinges for extra stability.

### Features

A slip-on knee brace.



# KNEE BRACES & SUPPORTS

Kenad SG Medical, Inc.

## Neoprene Knee Brace w/Anterior Closure



Size	Hinged Catalog No.	Stays Catalog No.
XSmall	16-3103	16-3203
Small	16-3104	16-3204
Medium	16-3105	16-3205
Large	16-3106	16-3206
XLarge	16-3107	16-3207
XXLarge	16-3108	16-3208
XXXLarge	16-3109	16-3209

### Construction

Royal Blue <sup>3</sup>/<sub>16</sub>" neoprene, nylon on both sides with hook and loop closure. Available with hinges or spiral stays.

### Features

Provides warmth, support, and equal compression. Hinges or spiral stays for extra stability. Popliteal cutout for patient comfort.

## Elastic Spiral Knee Brace



Size	Catalog No.
Small	16-4204
Medium	16-4205
Large	16-4206
XLarge	16-4207

### Construction

Medial and lateral spiral spring steel stays. Heavy two-ply elastic nylon/rubber blend.

### Features

Ideal for strains, sprains, and pre-op and post-op support.

## Patella Tendon Strap



Size	Catalog No.
Universal	16-6000

### Construction

Vinyl with neoprene sewn inside for greater compression of the patella tendon.

### Features

Contoured design allows full knee motion.

## Neoprene Knee Wrap w/Felt Buttress



Size	Catalog No.
Universal	16-8000

### Construction

<sup>1</sup>/<sub>8</sub>" neoprene hook compatible brushed nylon with universal felt buttress.

### Features

Provides warmth, support, and equal compression. Hook and loop closure for easy application.

# KNEE BRACES & SUPPORTS

Kenad SG Medical, Inc.

## Elastic Hinged Clinch Brace



Size	Catalog No.
Small	16-8004
Medium	16-8005
Large	16-8006
XLarge	16-8007

### Construction

High quality heavy two-ply elastic with loop/lock closure. Medial and lateral hinges provide rotational control and support.

## Neoprene Knee Wrap



Size	Catalog No.	Style
Adult	16-9000	Patella in
Child	16-9010	patella in
Adult	16-9100	Patella out
Child	16-9110	Patella out

### Construction

Neoprene rubber with nylon on both sides. Hook and loop closure.

### Features

Universal wrap-around design with popliteal opening for full knee flexion. Used in conservative treatment of arthritis, sprains, strains, or post-operative care. Available with patella in or out. Universal size reduces inventory requirements.

# ANKLE SUPPORTS

Kenad SG Medical, Inc.

## X-Brace Ankle Stirrup w/Foam Pad



**Size**  
Left  
Right

**Catalog No.**  
17-0020L  
17-0200R

### Construction

Hard plastic shell.

### Features

Figure "8" x-strap allows purchase on the instep and everts the foot. Also allows the strapping of the mid-foot to be tied to the stirrups.

## Gel/Air Ankle Stirrup



**Size**  
Small  
Regular  
Trainer

**Catalog No.**  
17-0100  
17-0110  
17-0120

### Construction

Hard plastic shell with air/gel bladder.

Hook and loop closure.

### Features

Combination air/gel bladder insert can be placed in freezer for two hours for cold therapy. Great for grades 1, 2, and 3 sprains when casting is not indicated. Fits right or left ankle to reduce inventory requirements.

## Short Leg Walker, STD



**Description**  
Small  
Medium/Large  
XLarge

**Catalog No.**  
17-0204  
17-0206  
17-0207

### Construction

Lightweight composite material used in malleable uprights for maximum strength.

### Features

Low profile rocker sole. Biomechanical design ensures normal walking gait. Heel design incorporates energy absorbing air-cells to reduce shock during heel strike.

## Deluxe Short Leg Walker



**Description**  
Small  
Medium  
Large

**Catalog No.**  
17-0204D  
17-0205D  
17-0206D

### Construction

Lightweight composite material with steel uprights for maximum strength.  $\frac{3}{8}$ " hook compatible foam laminated to wicking material to absorb perspiration. Hook and loop strapping ensures maximum stability.

### Features

Biomechanical design ensures normal walking gait. Heel design incorporates energy absorbing air-cells to reduce shock during heel strike. Ideal for nondisplaced fractures where casting is not indicated, minor strains, sprains, or after post cast removal.



# ANKLE SUPPORTS

Kenad SG Medical, Inc.

## Short Leg Walker – Short



Size	Catalog No.
Small	17-0304
Medium/Large	17-0306
XLarge	17-0307

### Construction

Lightweight composite material used in malleable uprights for maximum strength.

### Features

Low profile rocker sole. Biomechanical design ensures normal walking gait. Heel design incorporates energy absorbing air-cells to reduce shock during heel strike.

## Neoprene Ankle Sleeve



Size	Catalog No.
Small	17-0404
Medium	17-0405
Large	17-0406
XLarge	17-0407

### Construction

1/8" neoprene nylon, two sides.

### Features

Provides warmth, support, and equal compression.

## 90° Plantar Fasciitis Night Splint



Size	Catalog No.
Small	17-0504
Medium	17-0505
Large	17-0506

### Construction

Polyethylene shell with 1/8" foam laminated to tricot backing. Closed cell foam wedge. Hook and loop closure.

### Features

Tough flexible polyethylene shell cuts neatly and smoothly for customization. Air circulation and foot comfort are uncompromised. Wide strap pads help prevent strap pressure. Dorsiflexion wedge is included for additional stretching.

## Leather Lacing Ankle Support



Size	Catalog No.
1-2	17-1001
2-3	17-1002
3-4	17-1003
4-5	17-1004
7-10	17-1005
10-12	17-1006
12-14	17-1007
13-up	17-1008

### Construction

Soft leather.

### Features

This lace-up support provides excellent support for Grades I and II sprains. Fits left or right ankle. Available in eight sizes to fit exact patient needs.



# ANKLE SUPPORTS

Kenad SG Medical, Inc.

## Leather Ankle w/Elastic Wrap



Size	Catalog No.
Small	17-1004E
Medium	17-1005E
Large	17-1006E
XLarge	17-1007E

### Construction

Soft leather.

### Features

Lace-up support provides excellent support for Grade I and I ankle sprains. Elastic strap wraps in figure "8" for added stability.

## Leather Lacing Ankle Brace



Size	Catalog No.	Ankle Measurement
Small	17-1104	7.0" - 8.5"
Medium	17-1105	8.5" - 10.0"
Large	17-1106	10.0" - 11.5"
XLarge	17-1107	11.5" - 13.0"

### Construction

Laced leather support with removable aluminum stays.

### Features

Fits left or right ankle. Measure around smallest part of ankle (above malleolus) to determine correct size.

## Perforated Vinyl Ankle Brace w/Stays



Size	Catalog No.
Small	17-1204
Medium	17-1205
Large	17-1206
XLarge	17-1207

### Construction

Perforated vinyl laminate. Preformed aluminum stays.

### Features

Preformed ankle stays fits contour of the ankle. Fits either right or left ankle. Wipes clean with damp cloth.

## Nylon Ankle Stabilizer



Size	Catalog No.
Small	17-1304
Medium	17-1305
Large	17-1306
XLarge	17-1307

### Construction

Non-stretch nylon denier.

### Features

Eliminates taping. Hook and loop closure gives more adjustment capabilities. Available in black.

# ANKLE SUPPORTS

Kenad SG Medical, Inc.

## Vinyl Ankle Brace



Size	Catalog No.
Small	17-1504
Medium	17-1505
Large	17-1506
XLarge	17-1507

### Construction

Foam-vinyl laminate. Preformed aluminum stays fit contour of the ankle. Hook and loop closure.

### Features

Wipes clean. Fits left or right ankle.

## Canvas Lacing Ankle Brace



Size	Catalog No.	Ankle Measurement
Small	17-2004	7.0" - 8.5"
Medium	17-2005	8.5" - 10.0"
Large	17-2006	10.0" - 11.5"
XLarge	17-2007	11.5" - 13.0"

### Construction

Laced fabric support with removable aluminum stays.

### Features

Fits left or right ankle. Measure around smallest part of ankle (above malleolus) to determine correct size.

## Figure 8 Ankle Wrap



Size	Catalog No.
Small	17-2504
Medium	17-2505
Large	17-2506
XLarge	17-2507

### Construction

Heavy 3" elastic material.

### Features

Provides support for strains and sprains of the ankle.

## Tuli Heel Cups



Size	Catalog No.
Regular Rubber Under 175 lbs.	17-2604
Large Rubber Over 175 lbs.	17-2605
Regular Gel Under 175 lbs.	17-2606
Large Gel Over 175 lbs.	17-2607

### Construction

Molder (latex free) rubber.

### Features

New double-ribbed waffle design provides 50 percent more protection than previous model. Heel cups help prevent sore heels and knees, heel spurs, and shin splints. Available in two styles, molded rubber and polymer gel.

# ANKLE SUPPORTS

Kenad SG Medical, Inc.

## Visco-Elastic Heel Spur Pad



Men's Shoe Size	Women's Shoe Size	Catalog No.
2, 3, 4	4, 5, 6	17-2604A
5, 6, 7	7, 8, 9	17-2605B
8, 9, 10	10, 11, 12	17-2606C
11, 12, 13		17-2607D

### Construction

Made of visco-elastic polymer for optimum shock absorption.

### Features

Helps alleviate pain associated with heel spurs, also helps relieve the impact forces of heel strike. Designed to relieve pressure on muscles, tendons and joints in the lower extremities.

## Visco-Elastic Full Insoles



Men's Shoe Size	Women's Shoe Size	Catalog No.
2, 3, 4	4, 5, 6	17-2704A
5, 6, 7	7, 8, 9	17-2705B
8, 9, 10	10, 11, 12	17-2706C
11, 12, 13		17-2707D

### Construction

Visco-elastic polymer material for optimum shock absorption.

### Features

Help relieve the impact forces of heel strike. Also designed to relieve pressure on muscles, tendons, and joints in the lower extremities.

## Silicone Heel Cup w/Blue Spot



Men's Shoe Size	Women's Shoe Size	Catalog No.
2 to 7	4 to 8	17-2805
8 to 13	9 to 12	17-2806

### Construction

Silicone.

### Features

Blue dot is softer than surrounding material. Relieves pressure associated with heel spurs.

## Pull-On Ankle Support



Size	Catalog No.
Small	17-3604
Medium	17-3605
Large	17-3606
XLarge	17-3607

### Construction

Two-ply elastic.

### Features

Excellent for firm, gentle support.



# ANKLE SUPPORTS

Kenad SG Medical, Inc.

## Body Flex Ankle w/Figure "8" Strap



Size	Catalog No.
Small	17-3704
Medium	17-3705
Large	17-3706
XLarge	17-3707

### Construction

2.1mm thick nylon lycra.

### Features

Less heat, breathable, and thin as skin. Figure "8" strap provides greater stability greater compression, and non-slip suspension system.



# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Plaster Bandages & Splints



Description	Dimension	Catalog No.	Packaged
Elastic Plaster Rolls	4" x 5 Yds.	18-0104	Box of 12
Elastic Plaster Rolls	5" x 5 Yds.	18-0105	Box of 12
XF Plaster Rolls	2" x 3 Yds.	18-0202	Box of 12
XF Plaster Rolls	3" x 3 Yds.	18-0203	Box of 12
XF Plaster Rolls	4" x 5 Yds.	18-0204	Box of 12
XF Plaster Rolls	5" x 5 Yds.	18-0205	Box of 12
XF Plaster Rolls	6" x 5 Yds.	18-0206	Box of 12
XF Plaster Splints	3" x 15"	18-0303	Box of 50
XF Plaster Splints	4" x 15"	18-0304	Box of 50
XF Plaster Splints	5" x 30"	18-0305	Box of 50

### Features

Spooled rolls reduce telescoping. High quality plaster and substitute put more plaster in the cast and less on the floor. Available in elastic plaster rolls, extra fast rolls, and plaster splints.

## Plaster Prefab Splint Rolls



Description	Dimension	Catalog No.	Packaged
Maxi Splint Plaster Roll	2" x 15" x 24'	18-0402	Box of 2
Maxi Splint Plaster Roll	3" x 15" x 24'	18-0403	Box of 2
Maxi Splint Plaster Roll	4" x 15" x 24'	18-0404	Box of 2
Maxi Splint Plaster Roll	5" x 15" x 24'	18-0405	Box of 2
Maxi Splint Plaster Roll	6" x 15" x 24'	18-0406	Box of 2
Foam & Flannel Plaster Roll	2" x 15" x 20'	18-0502	Box of 2
Foam & Flannel Plaster Roll	3" x 15" x 20'	18-0503	Box of 2
Foam & Flannel Plaster Roll	4" x 15" x 20'	18-0504	Box of 2
Foam & Flannel Plaster Roll	5" x 15" x 20'	18-0505	Box of 2
Foam & Flannel Plaster Roll	6" x 15" x 20'	18-0506	Box of 2


### Construction

Maxi Splint is made of 1/4" polyurethane foam x-fast setting plaster of Paris encased in stockinett. Foam and flannel is made of higher density 3/8" green polyurethane form to ensure a more accurate fitting.

### Features

Available in 15 layers and is 24 feet long.

## Splint Cart Dispenser



Description	Catalog No.
Splint Cart Dispenser	18-1100

### Features

Designed to provide a lightweight alternative to bulky, heavy casting carts. Cart includes five compartment dispenser, waste basket and hanger, utility bucket and hanger. The perfect mobile casting and splinting cart for clinics that do not have a cast room.

# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Bandage Scissors



Size / Description	Catalog No.
Length 7 1/2" Double Ring	18-0600

### Construction

Steel chrome plated.

### Features

Surgical steel helps prevent rusting. Double ring for maximum torque when cutting multi-layers of material.

## Economy Vinyl Cast Boot



	Open Toe Catalog No.	Closed Toe Catalog No.	Dimension
Size			
XSmall	18-0703	18-0803	3" x 6 1/2"
Small	18-0704	18-0804	3 1/2" x 8 1/4"
Medium	18-0705	18-0805	4" x 9 1/2"
Large	18-0706	18-0806	4 1/4" x 11"
XLarge	18-0707	18-0807	4 1/2" x 12 1/2"

### Construction

Vinyl upper unit with hook and loop closures. Has 3/8" inner sole with cushioned outer sole.

### Features

Available in open- or closed-toe styles.

## Deluxe Canvas Rocker



Size	Catalog No.	Dimension
XXSmall	18-0902	3" x 5 1/5"
XSmall	18-0903	3 1/2" x 7 1/2"
Small	18-0904	4" x 9"
Medium	18-0905	4 1/4" x 10 1/2"
Large	18-0906	4 1/2" x 12"
XLarge	18-0907	5 5/8" x 13 1/2"

### Construction

Strong blue duck upper material. Hook and loop closures Cushioned inner Sole. No-skid rocker outer sole.

### Features

Rocker sole provides a more natural gait. May be worn over cast or as a shoe after foot surgery.

## Cast Toe Covers



Style	Catalog No.
Without Ankle Strap	18-1000
With Ankle Strap	18-1500

### Construction

Double layered stockinette. Available with elastic strap that slips on posteriorly or with elastic strap sewn into the end of the toe cover.

### Features

Ideal protection for toes in cold weather and to help keep toes clean. Universal size.



# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Economy Canvas Rocker



Size	Catalog No.
XXSmall	18-1002
XSmall	18-1003
Small	18-1004
Medium	18-1005
Large	18-1006
XLarge	18-1007

### Construction

Canvas upper with hook and loop closure. No-skid outer sole.

### Features

Rocker bottom provides a more natural gait.

## Deluxe Vinyl Cast Boot



Size	Open Toe Catalog No.	Closed Toe Catalog No.
XSmall	18-1203	18-1503
Small	18-1204	18-1504
Medium	18-1205	18-1505
Large	18-1206	18-1506
XLarge	18-1207	18-1507

### Construction

High quality upper vinyl unit. No-skid outer sole.

### Features

Rocker bottom provides a more natural gait. Available in open-toe and closed-toe styles.

## Shook Cast Stand



Description	Catalog No.
Shook Cast Stand	18-2000

### Construction

Stand is made of rust resistant chrome plated steel. Polyester foam pad and plastic bucket.

### Features

Excellent for maintaining ankle at 90 degrees during short leg casting.

## Rubber Walking Heel



Size	Catalog No.
Small	18-2004
Medium	18-2005
Large	18-2006

### Construction

An all-rubber walking heel with non-skid tread design.

### Features

A street-type design that is slightly contoured for easier walking. Easily incorporated in cast via ridges at each end to ensure firm seating. Sold in packages of 12.



# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Dual Walking Heel



**Size**  
Adult  
Child

**Catalog No.**  
18-2104  
18-2105

### Construction

Non-slip tread design.

### Features

Provides good patient stability by allowing regular heel-toe walking. Sold in dozen lots only.

## Post-Op Healing Shoe



### Men's Size

XSmall 18-2203  
Small 18-2204  
Medium 18-2205  
Large 18-2206  
XLarge 18-2207

### Catalog No.

### Women's Size

XSmall 18-2303  
Small 18-2304  
Medium 18-2305  
Large 18-2306  
XLarge 18-2307

### Catalog No.

### Construction

Breathable upper portion made of  $\frac{3}{8}$ " foam laminate. Hook and loop closure. Non-skid rubber sole.

### Features

Lightweight breathable upper with non-skid rubber bottom. Open flat for ease of foot insertion. Fits left or right foot to help reduce inventory requirements. Available in men and women's sizes.

## Cylinder Cast Support



### Style

Universal 18-2500  
LMC, Universal 18-2600

### Catalog No.

### Construction

Cotton webbing.

### Features

Quick release buckle for ease of removal. Prevents up and down pistoning motion of cylinder cast.

# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Deluxe Post-Op Shoe



Men's Size	Lacing Closure Catalog No.	Hook/Loop Closure Catalog No.
Small	18-2604	18-2504
Medium	18-2605	18-2505
Large	18-2606	18-2506
XLarge	18-2607	18-2507

Women's Size	Lacing Closure Catalog No.	Hook/Loop Closure Catalog No.
XSmall	18-3103	18-3003
Small	18-3104	18-3004
Medium	18-3105	18-3005
Large	18-3106	18-3006

### Construction

Nylon/foam upper unit, wood sole.

### Features

Available with lacing or hook and loop closures. Women's and men's styles.

## Economy Post-Op Shoe



Men's Size	Lacing Closure Catalog No.	Hook/Loop Closure Catalog No.	Dimension
Small	18-2704	18-2404	4" x 11 <sup>1</sup> / <sub>4</sub> "
Medium	18-2705	18-2405	4 <sup>1</sup> / <sub>4</sub> " x 11 <sup>1</sup> / <sub>2</sub> "
Large	18-2706	18-2406	4 <sup>1</sup> / <sub>2</sub> " x 12"
XLarge	18-2707	18-2407	4 <sup>1</sup> / <sub>2</sub> " x 12 <sup>3</sup> / <sub>4</sub> "

Women's Size	Lacing Closure Catalog No.	Hook/Loop Closure Catalog No.	Dimension
XSmall	18-2803	18-2903	3 <sup>1</sup> / <sub>4</sub> " x 9"
Small	18-2804	18-2904	3 <sup>1</sup> / <sub>2</sub> " x 9 <sup>1</sup> / <sub>2</sub> "
Medium	18-2805	18-2905	3 <sup>3</sup> / <sub>4</sub> " x 10"
Large	18-2806	18-2906	4" x 10 <sup>1</sup> / <sub>2</sub> "

### Construction

Strong vinyl upper unit, rigid wooden base with skid resistant ribbed rubber outer sole, with <sup>1</sup>/<sub>8</sub>" soft foam inner sole.

### Features

Available with lacing or hook and loop closures. Women's and men's styles.

# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## The Darco Post-Op Shoe



Men's Size	Catalog No.	Women's Size	Catalog No.
6 to 8	18-3204	4 to 6	18-3304
8 1/2 to 10	18-3205	6 1/2 to 8	18-3305
10 1/2 to 12	18-3206	8 1/2 to 10	18-3306
12 1/2 to 14	18-2206		
	18-3207		

### Construction

Soft padded mesh upper. Semi-ridged rubber sole. Hook and loop closure.

### Features

The Darco Shoe is generously proportioned to accommodate edema and foot dressings. Lightweight breathable upper with non-skid rubber bottom. Open flat for ease of foot insertion. Fits left or right foot to help reduce inventory requirements. Available in men's and women's sizes.

## Cotton Stockinette



Size	Catalog No.
2" x 25 yds.	18-3502
3" x 25 yds.	18-3503
4" x 25 yds.	18-3504

### Construction

Cotton.

### Features

Used to protect skin underneath plaster cast and also as a covering over a fiberglass cast. Can be used to hold bandage in place. Available in 2", 3", and 4" widths.

## Synthetic Stockinette



Size	Catalog No.
2" x 25 yds.	18-4002
3" x 25 yds.	18-4003
4" x 25 yds.	18-4004

### Construction

Polyester material.

### Features

Used to protect skin underneath plaster cast and also as a covering over a fiberglass cast. Can be used to hold bandage in place. Available in 2", 3", and 4" widths.



# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Orthopedic Felt



Size	Catalog No.
1/4" x 21" x 36"	18-4512
1/2" x 21" x 36"	18-4513

### Construction

Polyester fibers.

### Features

Will not retain water. Ideal for padding bony prominences under cast, neck pads, or arm slings, or padding of splint. Available in 1/2" and 1/4" thicknesses.

## Adhesive Tabbed Velfoam Straps



Size	Catalog No.
1" x 16"	18-5010
1" x 25"	18-5011
2" x 16"	18-5012
2" x 25"	18-5013

### Construction

3/8" velfoam with an adhesive pile tab at one end and a hook tab at the other end. Hook and loop closure.

### Features

Ideal for attaching to bivalved casts and splints to allow patient to remove and replace. Sold in dozen lots only. Box of 12.

## Adhesive Hook & Loop Rolls



Size	Type	Catalog No.
1" x 25 yds.	Loop	18-5307
1" x 25 yds.	Hook	18-5308
2" x 25 yds.	Loop	18-5309
2" x 25 yds.	Hook	18-5310

### Construction

Adhesive backed hook and loop rolls.

### Features

A must for cast rooms. Excellent for attaching to casts and splints. Available in 1" and 2" widths.

## Lewin Stern Splint



Type	Catalog No.
Finger	18-5407
Thumb	18-5408

### Construction

Available in padded style only.

### Features

Easily molded to finger, yet retains its strength. Used on fractures of the distal phalanx.

# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Hook & Loop Extension Tabs



**Catalog No.**  
18-5600

### Construction

Hook and loop extension tabs.

### Features

Ideal for extending strap when fitting obese patients with knee splints, slings, and wrist and forearm splints.

## Burnham Thumb & Finger Splines



### Type

Thumb	
Splint, Left	18-5804
Thumb	
Splint, Right	18-5805
Finger Splint	18-5806

### Catalog No.

### Construction

Plastic-coated aluminum. Padded.

### Features

Securely immobilize the wrist to arrest muscle-tendon tensions that tend to distort fracture fragments. Thumb splints come in left and right. Finger splint is pivotable.

## Elastic Bandages



### Size

2"	18-6010
3"	18-6011
4"	18-6012
6"	18-6013

### Catalog No.

### Construction

High quality rubber material.

### Features

Provides excellent support and compression to most any area of the body. Available in 2", 3", 4", and 6" widths. Sold in dozen lots only.

## Tubular Foam Padding



### Style

Blue Roll	18-6510
White Roll	18-6610

### Catalog No.

### Construction

3/8" foam with hook compatible brushed nylon laminate.

### Features

Fits over 1 1/2" or 1" webbing with buckles. Available in blue or white. Rolls are 12 feet long.

# CAST ROOM PRODUCTS

Kenad SG Medical, Inc.

## Bunion Splints



Size	Left/Right	Catalog No.
Small	Left	18-6704
Medium	Left	18-6705
Large	Left	18-6706
XLarge	Left	18-6707
Small	Right	18-6804
Medium	Right	18-6805
Large	Right	18-6806
XLarge	Right	18-6807

### Construction

Medial stay. Hook and loop closure.

### Features

Maintains great toe in neutral position at night while sleeping. Aluminum stay can be pre bent for custom toe alignment.

## Hanging Cast Sling



Size	Catalog No.
Small	18-7001
Medium	18-7002
Large	18-7003

### Features

This unique design features a quick release buckle system.



# SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Neoprene Tennis Elbow Sleeve w/Strap



Size	Catalog No.
XSmall	19-0403
Small	19-0404
Medium	19-0405
Large	19-0406
XLarge	19-0407

### Construction

1/8" neoprene nylon, two sides.

### Features

Provides warmth, support, and equal compression. Circumferential 2" hook and loop strap. Fits right or left forearm.

## Tennis Elbow Support w/Pad & Strap



Size	Catalog No.
XSmall	19-0503
Small	19-0504
Medium	19-0505
Large	19-0506
XLarge	19-0507
XXLarge	19-0508

### Construction

1/8" neoprene nylon, both sides. Elastic with loop lock closure.

### Features

4 1/2" circumferential neoprene sleeve tapering down to 3 1/2" width. Firm density foam pad provides greater compression to a smaller more focused area. Fits right or left arm. Ideal for the treatment of tennis elbow.

## Tennis Elbow Support



Description	Catalog No.
Link Closure	19-1000
Hook and Loop Closure	19-1001

### Construction

Made of 2" canvas with foam lining. Loop/lock or hook and loop closure.

### Features

Universal sizes.

## Body Flex Tennis Elbow Sleeve



Size	Catalog No.
Small	19-1104
Medium	19-1105
Large	19-1106
XLarge	19-1107

### Construction

2.1mm thick nylon lycra.

### Features

Less heat, breathable and thin as skin. Greater compression, and non-slip suspension system circumferential 2" hook and loop strap. Fits right or left forearm.

# SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Body Flex Tennis Elbow Strap w/Gel Pad



Size	Catalog No.
XSmall	19-1203
Small	19-1204
Medium	19-1205
Large	19-1206
XLarge	19-1207
XXLarge	19-1208

### Construction

2.1mm thick nylon lycra gel pad.

### Features

Less heat, breathable and thin as skin. Gel Pad and loop lock closure assures that direct pressure is applied in a smaller, focused area. Greater compression, and non-slip suspension system circumferential 2" hook and loop strap. Fits right or left forearm.

## Contoured Canvas Tennis Elbow



Size	Catalog No.
Small	19-2004
Medium	19-2005
Large	19-2006
XLarge	19-2007

### Construction

Made of 2" canvas with foam lining. Loop/lock closure.

### Features

Closure facilitates removal and replacement.

## Vinyl Tennis Elbow Splint



Size	Catalog No.
XSmall	19-2103
Small	19-2104
Medium	19-2105
Large	19-2106
XLarge	19-2107

### Construction

Made of 2" soft vinyl foam with loop/lock closure.

### Features

Closure permits patient to apply and remove the splint without assistance.

## Neoprene Pillow Tennis Elbow Support



Size	Catalog No.
Universal	19-3000

### Construction

$\frac{3}{8}$ " brushed pile velfoam with neoprene pillow incorporated into strap. Loop/lock closure.

### Features

Closure assures that direct pressure is applied in a smaller, focused area.

# SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Neoprene Tennis Elbow Strap



**Size**  
Universal

**Catalog No.**  
19-4000

**Construction**

1/8" neoprene hook compatible brushed nylon.

**Features**

Circumferential 2" hook and loop strap. Fits right or left forearm.

## Slip-on Elbow Support



**Size**  
XSmall  
Small  
Medium  
Large  
XLarge

**Catalog No.**  
19-4003  
19-4004  
19-4005  
19-4006  
19-4007

**Construction**

Two-ply elastic nylon/rubber blend material.

**Features**

Use where firm, gentle support is needed.

## Ulna Nerve Night Splint



**Size**  
Universal

**Catalog No.**  
19-6000

**Construction**

Blue foam with hook compatible pile. Hook closure.

**Features**

Designed to keep the elbow straight at night for patients that have nerve impingement when elbow is bent for long periods.



# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

Wrist splints and supports are used to provide temporary immobilization and support for stable fractures and ligamentous injuries to the wrist. They are also used to protect damaged soft tissue and for emergency immobilization prior to definitive treatment.

In emergency situations, the splint is often applied before the patient is moved to help reduce pain, improve circulation, and prevent further tissue injury. If there is an open fracture, there usually is no attempt to reposition, but a sterile dressing is applied and the arm is immobilized as gently as possible.

## 8" Cock Up Wrist Splint



Size	Catalog No.
Right Small	20-0204
Right Medium	20-0205
Right Large	20-0206
Left Small	20-0304
Left Medium	20-0305
Left Large	20-0306

### Construction

Brown vinyl laminated to  $\frac{3}{8}$ " foam with tricot lining. 1  $\frac{1}{2}$ " Volar stay.

### Features

Allows full motion of the MP's and limits wrist flexion. Loop lock closure assures easy adjustments.

## Plasti-Cast



Wrist Hand Size	Right Catalog No.	Left Catalog No.	Wrist Hand Thumb Size	Right Catalog No.	Left Catalog No.
Small	20-0404R	20-0404L	Small	20-0504R	20-0504L
Medium	20-0405R	20-0405L	Medium	20-0505R	20-0505L
Large	20-0406R	20-0406L	Large	20-0506R	20-0506L
XLarge	20-0407R	20-0407L	XLarge	20-0507R	20-0507L

### Construction

Semi-rigid polyethylene shell with a  $\frac{1}{4}$ " thermo-bonded closed cell foam. Loop lock closure.

### Features

Provides rigid wrist, hand, and thumb immobilization with full ROM. Ideal for minor strains, sprains, post operative, and post case removal.

## Vinyl Wrist & Forearm Splint



Size	Right 6" Catalog No.	Left 6" Catalog No.	Right 11" Catalog No.	Left 11" Catalog No.
XSmall	20-0803	20-0903	20-1103	20-1203
Small	20-0804	20-0904	20-1104	20-1204
Medium	20-0805	20-0905	20-1105	20-1205
Large	20-0806	20-0906	20-1106	20-1206
XLarge	20-0807	20-0907	20-1107	20-1207

### Construction

Brown vinyl laminated to  $\frac{3}{8}$ " foam with tricot lining.

### Features

Loop/lock closure at wrist provides greater support and compression.

# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Universal Wrist & Forearm Splint



Size	Catalog No.
XSmall	20-1003
Small	20-1004
Medium	20-1005
Large	20-1006

### Construction

Brown vinyl with  $\frac{3}{8}$ " foam backing.

### Features

Vinyl surface may be wiped clean. Universal design fits left or right hand.

## Short Perforated Vinyl Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XXSmall	20-1302	20-1402
XSmall	20-1303	20-1403
Small	20-1304	20-1404
Medium	20-1305	20-1405
Large	20-1306	20-1406
XLarge	20-1307	20-1407

### Construction

White perforated vinyl with soft cotton flannel backing.

### Features

Loop/lock closure at wrist for maximum wrist support. Cut out at palm allows for maximum carpal flexion.

## Contoured Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XXSmall	20-1502	20-1602
XSmall	20-1503	20-1603
Small	20-1504	20-1604
Medium	20-1505	20-1605
Large	20-1506	20-1606
XLarge	20-1507	20-1607

### Construction

White perforated vinyl with soft cotton flannel backing. Aluminum stays.

### Features

This splint is contoured at the wrist to ensure a more exact fit. Stays provide excellent support.

## Long Perforated Vinyl Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XXSmall	20-1702	20-1802
XSmall	20-1703	20-1803
Small	20-1704	20-1804
Medium	20-1705	20-1805
Large	20-1706	20-1806
XLarge	20-1707	20-1807

### Construction

White perforated vinyl with soft cotton flannel backing.

### Features

Loop/lock closure at wrist for maximum wrist support. Cut out at palm allows for maximum carpal flexion.



# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Contoured Canvas Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XXSmall	20-2102	20-2202
XSmall	20-2103	20-2203
Small	20-2104	20-2204
Medium	20-2105	20-2205
Large	20-2106	20-2206
XLarge	20-2107	20-2207

### Construction

Heavy blue canvas laminated to  $\frac{3}{8}$ " blue open-cell foam.

### Features

Loop/lock closure at wrist provides greater compression and support. Contoured at wrist to ensure a more accurate fit.

## Deluxe Universal Wrist & Forearm Splint



Length	Catalog No.
11"	20-2300

### Construction

$\frac{3}{8}$ " velfoam with hook and loop closures.

### Features

Dorsal stay pocket can be moved to allow a full range of size adjustments. Excess material may be trimmed for better fit. Universal design reduces inventory requirements.

## Trim Fit Universal Wrist & Forearm Splint



Length	Catalog No.
11"	20-2301 - Right 20-2302 - Left

### Construction

$\frac{3}{8}$ " velfoam with hook and loop closures.

### Features

Dorsal stay can be moved to allow a full range of size adjustments. Available in either right or left versions.

## Universal Elastic Wrist Support



Length	Catalog No.
XSmall	20-2403
Small	20-2404
Medium	20-2405
Large	20-2406
XLarge	20-2407

### Construction

Black elastic and two malleable aluminum stays.

### Features

Designed to fit either right or left hand by removing stays and turning over and reinserting stays. Universal design reduces inventory requirements.



# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Brown Vinyl Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XXSmall	20-2702	20-2802
XSmall	20-2703	20-2803
Small	20-2704	20-2804
Medium	20-2705	20-2805
Large	20-2706	20-2806
XLarge	20-2707	20-2807

### Construction

Brown vinyl laminated to  $\frac{3}{8}$ " foam with tricot lining.

### Features

Long wrist strap provides greater support. Wipes clean with damp cloth.

## Neoprene Wrist Splint w/Strap



Size	Catalog No.
Right Small	20-2904
Right Medium	20-2905
Right Large	20-2906
Right XLarge	20-2907
Left Small	20-3004
Left Medium	20-3005
Left Large	20-3006
Left XLarge	20-3007

### Construction

$\frac{1}{8}$ " neoprene nylon, two sides.

### Features

Provides warmth, support, and equal compression. Circumferential 3" hook and loop strap for added support with malleable volar stay.

## Blue Foam Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XSmall	20-3303	20-3403
Small	20-3304	20-3404
Medium	20-3305	20-3405
Large	20-3306	20-3406

### Construction

Brushed pile nylon laminated to blue foam.

### Features

Loop/lock closures. Long wrist strap provides greater support. Wipes clean with damp cloth.

# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Lacing Wrist Support



Size	Right 6" Catalog No.	Left 6" Catalog No.	Right 8" Catalog No.	Left 8" Catalog No.
XSmall	20-3703	20-3803	20-3903	20-4003
Small	20-3704	20-3804	20-3904	20-4004
Medium	20-3705	20-3805	20-3905	20-4005
Large	20-3706	20-3806	20-3906	20-4006
XLarge	20-3707	20-3807	20-3907	20-4007

### Construction

Hook/loop closure. Rigid spoon.

### Features

Closure permits application or removal with one hand. Once laces are tightened, hook and loop holds the splint firmly in place. Rigid spoon inside keeps hand in the dorsiflexed position. Spoon is removable and splint is washable.

## Pediatric Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
XXSmall	20-4502	20-4602
XSmall	20-4503	20-4603
Small	20-4504	20-4604

### Construction

Brushed pile nylon laminated to 3/8" blue foam.

### Features

Loop/lock closures. Long wrist strap provides greater support. Wipes clean with damp cloth. Available in three sizes to fit a wide range of pediatric patients.

## Neoprene Thumb & Wrist Strap



Description	Size	Catalog No.
Without Stay	Universal	20-4900
With Stay	Universal	20-5900

### Construction

Neoprene rubber with nylon on both sides.

### Features

Provides support and compression to the thumb and wrist. Available with or without stay.

# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Elastic Wrist & Forearm Splint



Size	Right Catalog No.	Left Catalog No.
Small	20-5104	20-5204
Medium	20-5105	20-5205
Large	20-5106	20-5206
XLarge	20-5107	20-5207

### Construction

High-quality two-ply elastic with nylon/cotton blend.

### Features

All-elastic design with hook closure allows firm support to the wrist.

## Vinyl Wrist & Forearm Splint w/Thumb Abducted



Size	Right Catalog No.	Left Catalog No.
Small	20-5304	20-5404
Medium	20-5305	20-5405
Large	20-5306	20-5406

### Construction

Brown vinyl laminated to  $\frac{3}{8}$ " foam with tricot lining.

### Features

Immobilizes thumb and wrist after strains, fractures, and dislocations.

## Canvas Wrist Brace w/ABD Thumb



Size	Right Catalog No.	Left Catalog No.
XSmall	20-5503	20-5603
Small	20-5504	20-5604
Medium	20-5505	20-5605
Large	20-5506	20-5606
XLarge	20-5507	20-5607

### Construction

$\frac{3}{8}$ " canvas foam laminate with cotton backing. Rigid thumb stay provides for immobilization of thumb. Loop lock closure.

### Features

Ideal for immobilization of the MP and CMC joint of the thumb. Allows for full MP flexion of the fingers. Volar stay adds greater stability at the wrist.

## Elastic Carpal Brace



Size	Catalog No.
Right XSmall	20-6203
Right Small	20-6204
Right Medium	20-6205
Right Large	20-6206
Right XLarge	20-6207
XSmall	20-6303
Left Small	20-6304
Left Medium	20-6305
Left Large	20-6306
XLarge	20-6307

### Construction

Perforated elastic with flexible spring steel stays.

### Features

Designed to help prevent carpal tunnel. Keeps wrist in neutral position. Allows free thumb movement and full MP motion.



# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Elastic Carpal Brace



Size	Catalog No.
Right XSmall	20-7003
Right Small	20-7004
Right Medium	20-7005
Right Large	20-7006
Right XLarge	20-7007
Left XSmall	20-7103
Left Small	20-7104
Left Medium	20-7105
Left Large	20-7106
Left XLarge	20-7107

### Construction

Heavy gauge elastic with 1 1/2" malleable volar stay.

### Features

Elastic provides compression and support. Volar stay for added stability. Hook and loop elastic straps wrap around and allows for increased or decreased compression at wrist.

## Exact Fit Deluxe Wrist Splint



Size	Catalog No.
Right XSmall	21-0503
Right Small	21-0504
Right Medium	21-0505
Right Large	21-0506
Right XLarge	21-0507
Left XSmall	21-0603
Left Small	21-0604
Left Medium	21-0605
Left Large	21-0606
Left XLarge	21-0607

### Construction

Nylon Denier, 2mm neoprene and 1" foam. 1 1/2" malleable Volar stay.

### Features

One of our most popular wrist support. Designed to provide heat, compression, and support at the wrist. Loop lock closure provides for exact fit.

## Exact Fit Deluxe Wrist Splint w/ABD Thumb



Size	Catalog No.
Right Small	21-0704
Right Medium	21-0705
Right Large	21-0706
Right XLarge	21-0707
Left Small	21-0804
Left Medium	21-0805
Left Large	21-0806
Left XLarge	21-0807

### Construction

Nylon Denier, 2mm neoprene and 1" foam. 1 1/2" malleable Volar stay. 1" ABD Thumb stay.

### Features

Immobilizes thumb and wrist after strains, fractures, and dislocations.

## Foam Wrist & Forearm Splint



Right Catalog No.	Left Catalog No.
21-1000	21-1100

### Construction

Brushed pile nylon laminated to blue 3/8" foam.

### Features

Two 3" elastic straps with hook tabs provide secure closure with compression. Universal size reduces inventory requirements.

# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Neoprene Wrist & Forearm Splint



**Size**  
Short 8"  
Long 10"

**Right  
Catalog No.**  
21-1400  
21-1600

**Left  
Catalog No.**  
21-1500  
21-1700

### Construction

Neoprene rubber with nylon on both sides.

### Features

Closure consists of two 3" elastic straps with hook tabs.

## Elastic Wrist Band



**Size**  
Universal

**Catalog No.**  
21-2400

### Construction

3" heavy elastic and brushed pile material.

### Features

Hook tabs ensure secure closure. Provides maximum support to the wrist while leaving the hand and fingers free.

## Black Universal Neoprene Wrist Wrap



**Size**  
Universal

**Catalog No.**  
21-2500

### Construction

1/8" neoprene brushed nylon.

### Features

Provides warmth, support, and equal compression to the wrist.

## Bilateral Elastic Wristlet



**Size**  
Small  
Medium  
Large  
XLarge

**Catalog No.**  
21-2804  
21-2805  
21-2806  
21-2807

### Construction

Heavy elastic with hook and loop closure.

### Features

May be used for either left or right hand simply by removing the stay and turning it over. Overall length 6".



# WRIST SPLINTS & SUPPORTS

Kenad SG Medical, Inc.

## Visco-Elastic Antivibration Glove



Size	Catalog No.
w/Thumb Small	21-4004
w/Thumb Medium	21-4005
w/Thumb Large	21-4006
w/Thumb XLarge	21-4007
w/o Thumb Small	21-4104
w/o Thumb Medium	21-4105
w/o Thumb Large	21-4106
w/o Thumb XLarge	21-4107

### Construction

Heavy leather, Visco-elastic polymer pad.

### Features

Designed to absorb shock in working environments that require the use of high impact tools for long periods. Helps prevent repetitive stress injuries of the wrist hand and shoulder. Visco-elastic polymer sewn in palm. Available with or without thumb.

## Thumb Holder



Hand Size	Size	Catalog No.
Right XSmall	2 <sup>1</sup> / <sub>4</sub> " to 2 <sup>3</sup> / <sub>4</sub> "	21-5003
Right Small	2 <sup>3</sup> / <sub>4</sub> " to 3 <sup>1</sup> / <sub>4</sub> "	21-5004
Right Medium	3 <sup>1</sup> / <sub>4</sub> " to 3 <sup>3</sup> / <sub>4</sub> "	21-5005
Right Large	3 <sup>3</sup> / <sub>4</sub> " to 4 <sup>1</sup> / <sub>2</sub> "	21-5006
Left XSmall	2 <sup>1</sup> / <sub>4</sub> " to 2 <sup>3</sup> / <sub>4</sub> "	21-5103
Left Small	2 <sup>3</sup> / <sub>4</sub> " to 3 <sup>1</sup> / <sub>4</sub> "	21-5104
Left Medium	3 <sup>1</sup> / <sub>4</sub> " to 3 <sup>3</sup> / <sub>4</sub> "	21-5105
Left Large	3 <sup>3</sup> / <sub>4</sub> " to 4 <sup>1</sup> / <sub>2</sub> "	21-5106

### Construction

Heat moldable thermal plastic.

### Features

Great for gamekeepers' thumb. Immobilizes the MP joint to protect injured collateral ligaments. Maintains thumb in opposition.

## Universal Carpal Brace



Size	Catalog No.
XSmall	21-6003
Small	21-6004
Medium	21-6005
Brace Large	21-6006
XLarge	21-6007

### Construction

Perforated elastic and four flexible spring steel stays.

### Features

Designed to fit either right or left hand by turning support over. Perforated elastic allows greater air flow. Universal design reduces inventory requirements.

## Carpal Brace w/Viscolas Pad



Size	Catalog No.
Right XSmall	21-6003V
Right Small	21-6004V
Right Medium	21-6005V
Right Large	21-6006V
Right XLarge	21-6007V
Left XSmall	21-6103V
Left Small	21-6104V
Left Medium	21-6105V
Left Large	21-6106V
Left XLarge	21-6107V

### Construction

Perforated elastic with Visco-elastic polymer pad.

### Features

Design for use with high impact air tools. Less bulky than the antivibration glove. Perforated elastic allows greater air flow. Full visco-elastic pad sewn on volar side



# KENAD RIB BELTS

Kenad SG Medical, Inc.

Rib belts provide compression for the thorax after rib fracture or surgery. Available in elastic or comfortable foam material.

## RIB BELT SIZE GUIDE

Size	Men's Chest Size	Women's Chest Size
XSmall	18" - 24"	18" - 24"
Small	24" - 30"	24" - 30"
Medium	30" - 36"	30" - 36"
Large	36" - 42"	36" - 42"
XLarge	42" - 48"	42" - 48"
XXLarge	48" - 52"	48" - 52"

## Comfort Universal Rib Belt



Chest Size	Men's Catalog No.	Women's Catalog No.
28" - 52"	22-1000	22-1100

### Construction

6" elastic and 1/8" foam laminated to a strong pile material.

### Features

Material permits controlled expansion and comfortable non-slip padding. Universal style reduces inventory requirements while allowing a full range of size adjustments.

## Universal Rib Belt



Chest Size	Men's 6" Catalog No.	Women's 6" Catalog No.	Men's 8" Catalog No.	Women's 8" Catalog No.
28" - 52"	22-1200	22-1300	22-1600	22-1700


### Construction

Elastic material. Hook and loop closure.

### Features

This universal design reduces inventory requirements while allowing a full range of size adjustments. Size adjustments are made possible by location of vertical pile strips positioned to correspond to standard rib belt sizes.

## Sized Elastic Rib Belt



Chest Size	Size	Men's Catalog No.	Women's Catalog No.
24" - 30"	Small	22-1404	22-1504
30" - 36"	Medium	22-1405	22-1505
36" - 42"	Large	22-1406	22-1506
42" - 48"	XLarge	22-1407	22-1507
48" - 52"	XXLarge	22-1408	22-1508

### Construction

6" elastic with an easy-to-fasten hook and loop closure.

### Features

Closure ensures support and comfort.

# KENAD ABDOMINAL BINDERS

Kenad SG Medical, Inc.

Abdominal binders are used to support abdominal muscles weakened by repeated surgery or multiple pregnancies. Binders are often used preoperatively for abdominal hernias and postoperatively for hysterectomies or other abdominal surgery. *ABDOMINAL BINDERS SHOULD NOT BE USED DURING PREGNANCY.*

## Universal Comfort Binder



**Hip Size**  
28" - 50"

**Catalog No.**  
23-1000

### Construction

10" wide elastic combined with  $\frac{3}{8}$ " foam laminated to brushed pile material.

### Features

Allows adjustment to all size ranges.

## Universal Abdominal Binder



**Hip Size**  
28" - 50"

**10" Width  
Catalog No.**  
23-1400

**12" Width  
Catalog No.**  
23-1500

### Construction

Made completely from elastic. Hook and loop closure.

### Features

Four pile strips are positioned vertically to correspond to standard abdominal binder sizes.

## Sized Abdominal Binder



**Hip Size**

**Size**

**10" Width  
Catalog No.**

**12" Width  
Catalog No.**

28" - 32"

Small

23-1704

23-1804

32" - 36"

Medium

23-1705

23-1805

36" - 40"

Large

23-1706

23-1806

40" - 45"

XLarge

23-1707

23-1807

45" - 50"

XXLarge

23-1708

23-1808

### Construction

Elastic with two form fitting stays. Hook and loop closure.

### Features

Anatomically correct custom form support to provide back support. Designed to be roll resistant for patient comfort.

# KENAD LUMBOSACRAL SUPPORTS

Kenad SG Medical, Inc.

Lumbosacral supports help to relax the sacro-spinal muscles to modify pelvic inclination and lumbar lordosis, and provide support for lower back injuries. These supports are often used in conservative treatment of back problems and for postoperative lumbosacral support and post-op work hardening therapy programs. Many industries require employees who perform lifting tasks to wear special lumbosacral supports.

## Vinyl Lumbosacral Support



Size	11" Width Catalog No.	13" Width Catalog No.
XSmall	24-1003	24-1503
Small	24-1004	24-1504
Medium	24-1005	24-1505
Large	24-1006	24-1506
XLarge	24-1007	24-1507
XXLarge	24-1008	24-1508

### Construction

Quality elastic with  $\frac{3}{8}$ " foam laminated to vinyl in the back. Four contoured spring steel stays.

### Features

Stays provide extra rigid support. Hook and loop closure assures easy application and removal. Retaining straps provide additional support in front.

## Lumbosacral Support w/Foam Pad



Size	Catalog No.
Small	24-1204
Medium	24-1205
Large	24-1206
XLarge	24-1207
XXLarge	24-1208
XXXLarge	24-1209

### Construction

10" elastic with firm pad.

### Features

Elastic construction with pad in sacral area for added support. 3" side pulls allows maximum sacral compression.

## Lumbosacral Support



Hip Size	11" Width Size	Catalog No.
28" - 32"	XSmall	24-1103
32" - 36"	Small	24-1104
36" - 40"	Medium	24-1105
40" - 44"	Large	24-1106
44" - 48"	XLarge	24-1107

### Construction

Quality elastic with laminated  $\frac{3}{8}$ " foam pad in the back. Contoured stays. Hook and loop closure.

### Features

An extremely comfortable belt. Retaining straps provide additional support in front.



# KENAD LUMBOSACRAL SUPPORTS

Kenad SG Medical, Inc.

## Mesh Lumbosacral Support



Hip Size	11" Width Size	Catalog No.
28" - 48"	Universal	24-1400
28" - 32"	XSmall	24-1403
32" - 36"	Small	24-1404
36" - 40"	Medium	24-1405
40" - 44"	Large	24-1406
44" - 48"	XLarge	24-1407

### Construction

Open mesh back. Elastic with four metal stays in back. Hook and loop closure.

### Features

Open mesh back is well suited to wear in warm weather. Retaining straps provide additional support in front.

## Criss/Cross Lumbosacral Support



Hip Size	Size	Catalog No.
28" - 48"	Universal	24-1600
28" - 32"	XSmall	24-1603
32" - 36"	Small	24-1604
36" - 40"	Medium	24-1605
40" - 44"	Large	24-1606
44" - 48"	XLarge	24-1607
48" - 52"	XXLarge	24-1608

### Construction

All elastic construction features a criss/cross design in back with adjustable side pull straps to allow proper tension in the sacral area.

### Features

Unique support is 10" wide in back with hook and loop closure for easy application and removal.

## Lumbosacral w/Neoprene Pocket



Size	Catalog No.
XSmall	24-1703
Small	24-1704
Medium	24-1705
Large	24-1706
XLarge	24-1707
XXLarge	24-1708

### Construction

9" and 4" elastic, 2mm neoprene nylon both sides.

### Features

Elastic construction with neoprene pocket for thermal plastic insert for additional support. Four-inch sides pulls allows for greater compression at the sacral area.

## 10" Neoprene Lumbosacral w/Pocket



Size	Catalog No.
XSmall	24-1803
Small	24-1804
Medium	24-1805
Large	24-1806
XLarge	24-1807
XXLarge	24-1808

### Construction

1/8" neoprene nylon both sides.

### Features

All neoprene construction with neoprene pocket for thermal plastic insert for additional support. 4 1/2" sides pulls allows for greater compression at the sacral area. Retains heat to keep muscles relaxed.

# KENAD LUMBOSACRAL SUPPORTS

Kenad SG Medical, Inc.

## Body Flex 9" Lumbosacral Support w/Pocket



Size	Catalog No.
XSmall	24-1903
Small	24-1904
Medium	24-1905
Large	24-1906
XLarge	24-1907
XXLarge	24-1908

### Construction

2.1mm thick nylon Lycra and perforated elastic with Lycra pocket for thermal plastic insert. Four-inch side pulls allows for greater compression at the sacral area.

### Features

Less heat, breathable, and thin as skin. Retains heat to keep muscles relaxed.

## Chair Back Sling



Size	Catalog No.
Universal	24-2000

### Construction

Durable Tricot fabric with memory foam pad.

### Features

Designed to provide relief for low back pain. Prevents muscle fatigue and discomfort. Fits most chairs. Great for the home or office. Adjust to fit your lumbar curve. Removable memory foam pad may be placed in any position.

## Lumbar Roll



Description	Size	Catalog No.
Lumbar Roll	Universal	24-2300
Half Lumbar Roll		24-2350

### Construction

Medium density foam covered with a cotton/polyester cover.

### Features

Helps reduce back fatigue when sitting for long periods.

## Lumbosacral Cushion



Description	Catalog No.
Lumbosacral Cushion	24-2400

### Construction

Medium density foam covered with a cotton/polyester cover.

### Features

Provides support for mid-thoracic to sacrum area.

# KENAD LUMBOSACRAL SUPPORTS

Kenad SG Medical, Inc.

## Economy Lumbosacral Support



Size	Catalog No.
XSmall	24-2903
Small	24-2904
Medium	24-2905
Large	24-2906
XLarge	24-2907

### Construction

All elastic construction with adjustable side pull straps. Four spring steel stays.

### Features

Pull straps allow proper tension adjustment in the lumbar area. Stays provide maximum support. Sizes fit most patient needs.

## Lo-Bak™ Industrial Lumbosacral Support



Size	Straps Catalog No.	Detach Straps Catalog No.	No Straps Catalog No.
Small	24-3104V	24-4104V	24-5104V
Medium	24-3105V	24-4105V	24-5105V
Large	24-3106V	24-4106V	24-5106V
XLarge	24-3107V	24-4107V	24-5107V
XXLarge	24-3108V	24-4108V	24-5108V

### Construction

Sturdy elastic panels provide abdominal lift that helps reduce lower back fatigue. Black elastic material. Hook and loop closure. Double pull. Four strategically positioned stays. Available with or without permanent or removable adjustable shoulder straps.

### Features

Provides firm support to the lumbosacral portion of the spine for workers whose jobs involve lifting and lowering. Also useful in work hardening therapy postoperatively or after conservative treatment.

## Padded Thermoplastic Lumbar Insert



Size	Catalog No.
Universal	24-6000

### Construction

Thermoplastic material.

### Features

Designed for use in lumbosacral supports with pocket. Can be heated and molded for additional support. Pre-cut to conform to anatomical contours.



# KENAD LUMBOSACRAL SUPPORTS

Kenad SG Medical, Inc.

## Lumbosacral Cold Wrap



**Size**  
Universal

**Catalog No.**  
24-7000

### Construction

$\frac{3}{8}$ " brushed pile foam hook compatible.

### Features

Designed for use with cold therapy for the lumbar spine.

## Lo-Bak 10" Industrial Back Support



**Size**  
Small  
Medium  
Large  
XLarge  
XXLarge

**Catalog No.**  
24-7004  
24-7005  
24-7006  
24-7007  
24-7008

### Construction

Perforated elastic with loop and  $\frac{1}{8}$ " foam laminate.

### Features

Provides support to the lumbosacral portion of the spine for workers whose jobs involve lifting and lowering. Sturdy panels provide abdominal lift that helps reduce low back fatigue. Side pulls allow for added support. Comes with shoulder straps.

## Shoulder Cold Wrap



**Size**  
Universal

**Catalog No.**  
24-7010

### Construction

$\frac{3}{8}$ " brushed pile foam hook compatible.

### Features

Designed for use with cold therapy for the shoulder.

## Knee Cold Wrap



**Size**  
Universal

**Catalog No.**  
24-7020

### Construction

$\frac{3}{8}$ " brushed pile foam hook compatible.

### Features

Designed for use with cold therapy for the knee.

# KENAD LUMBOSACRAL SUPPORTS

Kenad SG Medical, Inc.

## Economy 6" Lo-Bak Industrial Back Support



Description	Size	Catalog No.
N/Strap Black	Small	24-8004
N/Strap Black	Medium	24-8005
N/Strap Black	Large	24-8006
N/Strap Black	XLarge	24-8007
Perm/Strap Black	Small	24-8104
Perm/Strap Black	Medium	24-8105
Perm/Strap Black	Large	24-8106
Perm/Strap Black	XLarge	24-8107
D/Strap Black	Small	24-8204
D/Strap Black	Medium	24-8205
D/Strap Black	Large	24-8206
D/Strap Black	XLarge	24-8207

### Construction

All elastic construction.

### Features

Provides support to the lumbosacral portion of the spine. For workers whose jobs involve lifting and lowering. Sturdy panels provide abdominal lift that helps reduce low back fatigue. Side pulls allow for added support. Comes with optional shoulder straps. Permanent, detachable, or no straps.

## Lycra Lo-Bak Industrial Back Support



Size	Catalog No.
XSmall	24-7803
Small	24-7804
Medium	24-7805
Large	24-7806
XLarge	24-7807
XXLarge	24-7808

### Construction

All mesh stretchable lycra and elastic.

### Features

Provides support to the lumbosacral portion of the spine for workers whose jobs involve lifting and lowering. Sturdy panels provide abdominal lift which helps reduce low back fatigue. Side pulls allow for added support. Comes with shoulder straps.

# CLAVICULAR SUPPORTS

Kenad SG Medical, Inc.

Used for clavicular fractures where there is no articular involvement. Also good for stooped shoulders.

## Foam Clavicle Brace



Size	Catalog No.
XXXSmall	25-1001
XXSmall	25-1002
XSmall	25-1003
Small	25-1004
Medium	25-1005
Large	25-1006
XLarge	25-1007
XXLarge	25-1008

### Construction

Cotton webbing covered with tube foam laminated to a brushed pile material that is hook compatible.

### Features

Hook and loop closure.

## Clavicle Brace



Size	Stretch Style Catalog No.	Firm Style Catalog No.
XSmall	25-1103	25-1203
Small	25-1104	25-1204
Medium	25-1105	25-1205
Large	25-1106	25-1206
XLarge	25-1107	25-1207

### Construction

Pads are made of 1/2" foam covered with stockinette.

### Features

Straps fasten to triangular buckle closure in back. Pronged buckles secure the brace.

## Heavy Duty Foam Clavicle



Catalog No.
25-1303
25-1304
25-1305
25-1306
25-1307

### Construction

Cotton stockinette and medium density foam. Hook and loop closure.

### Features

Designed for stabilization of clavicle fractures. Hook and loop closure for easy application and removal.



# CLAVICULAR SUPPORTS

Kenad SG Medical, Inc.

## Clavice Strap



### Description

Unpadded Clavice Strap  
Unpadded Clavice Strap  
Unpadded Clavice Strap  
Unpadded Clavice Strap  
Padded Clavice Strap  
Padded Clavice Strap  
Padded Clavice Strap  
Padded Clavice Strap

### Size

Small  
Medium  
Large  
XLarge  
Small  
Medium  
Large  
XLarge

### Catalog No.

25-1604  
25-1605  
25-1606  
25-1607  
25-1704  
25-1705  
25-1706  
25-1707

### Construction

1 1/2" poly/cotton webbing with plastic retainer for straps.

### Features

Designed for stabilization of clavicle fractures. Strap closes in front with metal buckle for easy application and removal.

## Padded Clavicle & Shoulder Brace



### Size

XSmall  
Small  
Medium  
Large  
XLarge

### Catalog No.

25-1803  
25-1804  
25-1805  
25-1806  
25-1807

### Construction

Axilla pads are medium density foam covered with stockinette.

### Features

An exceptionally well padded brace of one-piece construction. Designed to assure positive support and patient comfort. Adjustment straps hold axilla pads up and back, providing good control of forward and downward displacement of the shoulder. Slide buckle and catch provide one-time adjustment and facilitates removal.

## Padded Clavicle Brace



### Size

XSmall  
Small  
Medium  
Large  
XLarge

### Catalog No.

25-1903  
25-1904  
25-1905  
25-1906  
25-1907

### Construction

Axilla pads are 3/8" foam laminate.

### Features

Well padded for conditions involving scapular and clavicular areas. Each strap is fastened on both ends by a two-prong buckle to provide a snug, no-slip fit. Four adjustment points permit a more exact fit.

# CLAVICULAR SUPPORTS

Kenad SG Medical, Inc.

## Heavy Duty Felt Clavicle



Size	Catalog No.
XSmall	25-2003
Small	25-2004
Medium	25-2005
Large	25-2006
XLarge	25-2007

### Construction

2" webbing, orthopedic felt, and stockinette.

### Features

Designed for stabilization of clavicle fractures. Stockinette covered felt. Felt padding helps prevent pressure points.

# CLAVICULAR SUPPORTS

Kenad SG Medical, Inc.

## ACROMIOCLAVICULAR BRACES

Used in the treatment of acromioclavicular separations and dislocations. Provide steady upward pressure on the arm from the elbow to the dislocation site and downward pressure on the acromioclavicular region.

### Hard Shell Acromioclavicular Brace



Size	Catalog No.
Small	25-3504
Large	25-3506

#### Construction

$\frac{3}{8}$ " foam laminated to vinyl and covered with royal blue cotton duck material attached to an aluminum posterior elbow splint.

#### Features

Quick release buckle system facilitates application and removal.

### Acromioclavicular Brace – Short



Size	Catalog No.
Small	25-4004
Medium	25-4005
Large	25-4006

#### Construction

$\frac{3}{8}$ " foam laminated to vinyl and covered with royal blue cotton duck material. Soft shell. Short pocket.

#### Features

Quick release buckle system facilitates application and removal. Our most popular design.

### Acromioclavicular Brace – Long



Size	Catalog No.
Small	25-4504
Medium	25-4505
Large	25-4506

#### Construction

$\frac{3}{8}$ " foam laminated to vinyl and covered with royal blue cotton duck material. Soft shell. Long pocket.

#### Features

Longer pocket supports hand and wrist. Quick release buckle system allows easy application and removal.



# CLAVICULAR SUPPORTS

Kenad SG Medical, Inc.

## Universal Economy Acromioclavicular Brace

**Size**

Universal

**Catalog No.**

25-4600

**Construction**

Cotton duck material lined with strips of foam.  
Shoulder pad is vinyl laminated to  $\frac{3}{8}$ " foam. Tricot lined.

# CERVICAL COLLARS

Kenad SG Medical, Inc.

Available in straight, serpentine, and split collar styles and in universal and sized versions. Offered in four heights from 2" to 4 1/2" in soft, medium, and firm density foam covered with cotton stockinette.

## Cerve-Ease Collar



Dimension	Size	Catalog No.
3" x 19"	XSmall - XLong	26-2103
3" x 17"	Small	26-2104
4 1/2" x 19"	Medium	26-2105
4 1/2" x 21"	Large	26-2106

### Construction

Flexible plastic stiffener fastened to 1" foam. Hook and loop closure.

### Features

Comfortable, rigid support.

## Universal Collars



Medium Density Catalog No.	Firm Density Catalog No.	Dimension
26-2201	26-2301	2" x 23"
26-2202	26-2302	2 1/2" x 23"
26-2203	26-2303	3" x 23"
26-2204	26-2304	3 1/2" x 23"
26-2205	26-2305	4" x 23"
26-2206	26-2306	4 1/2" x 23"

### Construction

1" foam covered with stockinette. Medium and firm density foam. Hook and loop closure.

### Features

Designed for minor sprains, strains, or after post cast, or brace removal of the cervical spine.

## Split Cervical Collar



Size	Catalog No.	Dimension
Small	26-2601	4" x 26"
Medium	26-2602	4" x 28"
Large	26-2603	4" x 30"

### Construction

1" medium density foam covered with stockinette. Metal grommets used with a spreader bar attachment when used as a head halter. Hook and loop closure.

### Features

Versatile collar designed so that it can be used as a head halter.

## Philadelphia Cervical Collar®\* Liner



Description	Catalog No.
Philadelphia Cervical Collar* Liner	26-5000

### Construction

Two layers of cotton stockinette.

### Features

Used to line collar for patient skin protection.

*\*The Philadelphia Cervical Collar® is a registered trademark of Charles Greiner and Company, Inc.*

# CERVICAL COLLARS

Kenad SG Medical, Inc.

## Serpentine Collar



Collar Height	Length	Size	Medium Catalog No.	Firm Catalog No.
2"	12 <sup>1</sup> / <sub>2</sub> "	XS-Pedi	26-2401	26-2501
2"	21"	Medium	26-2402	26-2502
2"	24"	Large	26-2403	26-2503
2 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	Small	26-2404	26-2504
2 <sup>1</sup> / <sub>2</sub> "	18 <sup>1</sup> / <sub>2</sub> "	Medium	26-2405	26-2505
2 <sup>1</sup> / <sub>2</sub> "	20 <sup>1</sup> / <sub>2</sub> "	Large	26-2406	26-2506
2 <sup>1</sup> / <sub>2</sub> "	21 <sup>1</sup> / <sub>2</sub> "	XLarge	26-2407	26-2507
3"	15 <sup>1</sup> / <sub>2</sub> "	Small	26-2408	26-2508
3"	17"	Medium	26-2409	26-2509
3"	20"	Large	26-2410	26-2510
3"	21 <sup>1</sup> / <sub>2</sub> "	XLarge	26-2411	26-2511
3 <sup>3</sup> / <sub>4</sub> "	19"	Medium	26-2412	26-2512
3 <sup>3</sup> / <sub>4</sub> "	22"	XLong Med	26-2413	26-2513
4"	20 <sup>1</sup> / <sub>2</sub> "	Small	26-2414	26-2514
4"	21 <sup>1</sup> / <sub>2</sub> "	Medium	26-2415	26-2515
4"	24 <sup>1</sup> / <sub>2</sub> "	Large	26-2416	26-2516
4 <sup>1</sup> / <sub>2</sub> "	22"	Large	26-2417	26-2517
5"	21 <sup>1</sup> / <sub>2</sub> "	Large	26-2418	26-2518
5"	24"	XLarge	26-2419	26-2519

### Construction

Constructed from either medium or firm density 1" foam covered by cotton stockinette.

### Features

Designed to provide comfortable, semi-rigid immobilization. Shaped to conform to patient's chin or jaw. Extra stockinette cover included.

## Contoured Collar



Collar Height	Length	Size	Catalog No.
3"	16"	Small	26-2700
3"	18"	Medium	26-2701
3"	20"	Large	26-2702
3"	22"	XLarge	26-2703
3 <sup>1</sup> / <sub>2</sub> "	16"	Small	26-2704
3 <sup>1</sup> / <sub>2</sub> "	18"	Medium	26-2705
3 <sup>1</sup> / <sub>2</sub> "	20"	Large	26-2706
3 <sup>1</sup> / <sub>2</sub> "	22"	XLarge	26-2707
4"	16"	Small	26-2708
4"	18"	Medium	26-2709
4"	20"	Large	26-2710
4"	22"	XLarge	26-2711

### Construction

Medium density 1" foam covered with stockinette.

### Features

Hook and loop closure.



# CERVICAL COLLARS

Kenad SG Medical, Inc.

## Philadelphia Cervical Collar®\*



Collar Height	Size	Catalog No.
2 1/4"	Small	26-5004
2 1/4"	Medium	26-5005
2 1/4"	Large	26-5006
3 1/4"	Small	26-5104
3 1/4"	Medium	26-5105
3 1/4"	Large	26-5106
4 1/4"	Small	26-5204
4 1/4"	Medium	26-5205
4 1/4"	Large	26-5206
5 1/4"	Small	26-5304
5 1/4"	Medium	26-5305
5 1/4"	Large	26-5306

### Construction

Plastazote foam with hook and loop fasteners.

### Features

The original Philadelphia Cervical Collar®. Used to limit cervical motion and help reduce nerve root irritation. Occipital and mandibular posts provide stability. May be utilized in cervical casting procedures as well as cervical immobilization in emergency situations. Four collar heights provide flexibility in achieving flexion, neutral position, or hyperextension.

*\*The Philadelphia Cervical Collar® is a registered trademark of Charles Greiner and Company, Inc.*

## Exo-Static Collar



Size	Catalog No.
Small	26-6104
Small-Long	26-6104SL
Medium	26-6105
Large	26-6106

### Construction

Lightweight flexible plastic.

### Features

One-piece design allows for easy application. Height is adjustable by using the telescoping slotted metal bars.

## Exo-Static Collar w/Chin Piece



Size	Catalog No.
Small	26-6204
Medium	26-6205
Large	26-6206

### Construction

Lightweight flexible plastic.

### Features

One-piece design allows for easy application. Height is adjustable by using the telescoping slotted metal bars. Chin piece adds greater stability.

## Philly A-Series Collar



Size	Catalog No.
Infant	26-7002I
Pediatric	26-7003P
Stout	26-7004S
Medium	26-7005M
Std.	26-7006STD
Tall	26-7007T

### Construction

Lightweight flexible plastic. Hook compatible padded foam.

### Features

Unique design reduces rotation and slightly hyperextends the cervical spine. Two-piece design provides ridged immobilization. Padded foam for added comfort. Hook and loop closure facilitates easy application and removable.

# KENAD SHOULDER IMMOBILIZERS, SLINGS

Kenad SG Medical, Inc.

Shoulder immobilizers have a wide range of uses from strains, subluxations and dislocations to rupture of the shoulder tendon, fractures of the scapula, clavicle and humerus and following shoulder joint replacement. All shoulder immobilizers have foam padded shoulder straps, unless otherwise noted.

## Deluxe Sized Shoulder Immobilizer



Size	Catalog No.
Small-Medium	27-1404
Large-XLarge	27-1406

### Construction

$\frac{3}{8}$ " foam laminated hook and loop closure.

### Features

Our most popular shoulder immobilizer. Humeral cuff, wrist cuff, and shoulder strap can be positioned where desired. Shoulder strap with buckle for easy adjustment.

## Economy Universal Shoulder Immobilizer



Size	Catalog No.
Adult	27-1500
Child	27-1510

### Construction

$\frac{3}{8}$ " foam laminated to a brushed pile nylon. Hook and loop closure.

### Features

Immobilize either arm, shoulder, or both during convalescence from dislocation, fracture, or surgery.

## Deluxe Universal Shoulder Immobilizer



Size	Catalog No.
Universal	27-1600

### Construction

$\frac{3}{8}$ " foam laminated to a brushed pile nylon. Padded wrist cuff and heavy non-stretch arm strap.

### Features

Hook and loop closure facilitates application and removal.

## Elastic Shoulder Immobilizer



Size	Men's Catalog No.	Women's Catalog No.
XSmall	27-2003	27-2503
Small	27-2004	27-2504
Medium	27-2005	27-2505
Large	27-2006	27-2506
XLarge	27-2007	27-2507
XXLarge	27-2008	27-2508

### Construction

Heavy elastic.

### Features

Hook and loop closure allows immobilization of shoulder and arm.



# KENAD SHOULDER IMMOBILIZERS, SLINGS

Kenad SG Medical, Inc.

## Deluxe Elastic Shoulder Immobilizer



### Adjustment

#### Range

24" - 30"

30" - 36"

36" - 42"

42" - 48"

#### Size

Small

Medium

Large

XLarge

#### Catalog No.

27-2704

27-2705

27-2706

27-2707

### Construction

Kodel and elastic with hook and loop closure.

### Features

Immobilizes either arm, shoulder, or both during convalescence from dislocation, fracture, or surgery.

## Arm Sling Immobilizer



### Size

XSmall

Small

Medium

Large

XLarge

### Catalog No.

27-2903

27-2904

27-2905

27-2906

27-2907

### Construction

Constructed from blue material. 1" webbing. Hook and D-ring closure.

### Features

Can be used effectively as either an arm sling or shoulder immobilizer.

## Velpeau Sling



### Size

XSmall

Small

Medium

Large

XLarge

### Catalog No.

27-3003

27-3004

27-3005

27-3006

27-3007

### Construction

Beige cotton duck material with cotton webbing.

### Features

Can be used as either an arm sling or shoulder immobilizer.

## Universal Sling & Swath



### Size

Adult

Child

### Catalog No.

27-3500

27-3503

### Construction

$\frac{3}{8}$ " velfoam.

### Features

Hook and loop closures permit a wide range of adjustment. Universal size reduces inventory requirements. Fits right or left arm.



# KENAD SHOULDER IMMOBILIZERS, SLINGS

Kenad SG Medical, Inc.

## Sling & Swath



### Adjustment

Range	Size	Catalog No.
24" x 30"	Small	27-3704
30" x 36"	Medium	27-3705
36" x 42"	Large	27-3706
42" x 48"	XLarge	27-3707

### Construction

Sling and swath wrap are made from  $\frac{3}{8}$ " velfoam.

### Features

Hook and loop closures permit a wide range of adjustment. Designed to fit either right or left arm.

## Humeral Abduction Pillow



### Size

Small  
Large

### Catalog No.

27-4004  
27-4006

### Construction

Made of firm density foam covered with  $\frac{3}{8}$ " velfoam.

### Features

Used in the treatment of shoulder injuries. Pillow can be turned to maintain 75-degree or 45-degree angle. Also comes with an extension pillow for a larger arm rest.

## Arm Sling with Swath



### Size

Adult

### Catalog No.

27-5000

### Construction

Made of blue poly/cotton material. Swath is made of  $\frac{3}{8}$ " velfoam with hook and loop closure.

### Features

Can be used as an arm sling or shoulder immobilizer. Universal size reduces inventory requirements.

# KENAD ARM SLINGS

Kenad SG Medical, Inc.

Used in the treatment of bone, muscular, or ligamentous injuries to the hand, arm, elbow, or shoulder. Also used to provide arm support after surgical or nonsurgical correction, especially in conjunction with a cast. Often used in the field or emergency room to immobilize a fracture or dislocation awaiting physician intervention. All slings have foam neck pads unless otherwise noted.

## Universal Arm Sling



Size	Beige Catalog No.	Blue Catalog No.
Adult	28-1400	28-1500
Child	28-1700	28-1800

## Economy Arm Sling



Size	Beige Catalog No.	Blue Catalog No.
XSmall	28-2003	28-8003
Small	28-2004	28-8004
Medium	28-2005	28-8005
Large	28-2006	28-8006
XLarge	28-2007	28-8007

### Construction

Beige or royal blue cotton duck material. 1" cotton webbing.

### Features

Sold in dozen lots only. No neck pad.

## Deluxe Arm Sling w/Foam Strap



Size	Catalog No.
XSmall	28-2203
Small	28-2204
Medium	28-2205
Large	28-2206
XLarge	28-2207

### Construction

$\frac{3}{8}$ " foam laminated to a brushed pile nylon. Padded wrist cuff and heavy non-stretch arm strap.

### Features

Hook and loop closure facilitates application and removal.

## Arm Sling



Size	Catalog No.
XXSmall	28-2502
XSmall	28-2503
Small	28-2504
Medium	28-2505
Large	28-2506
XLarge	28-2507

### Construction

Beige. 1 $\frac{1}{2}$ " cotton webbing. Hook and loop closure.

# KENAD ARM SLINGS

Kenad SG Medical, Inc.

## Nylon Arm Sling



Size	Catalog No.
Small	28-2604
Medium	28-2605
Large	28-2606
XLarge	28-2607

### Construction

Heavy black denier nylon. 1 1/2" heavy poly webbing.

### Features

Foam neck pad for patient comfort. Slide buckle on strap for easy adjustment. Thumb loop allows support of hand. High quality heavy material for durability and patient comfort. Hook and Loop closure.

## Denim Arm Sling



Size	Catalog No.
XSmall	28-3003
Small	28-3004
Medium	28-3005
Large	28-3006
XLarge	28-3007

### Construction

High quality heavy denim.

### Features

Wider shoulder strap helps distribute arm weight more evenly. Hook and loop closure with side buckle permits a wide range of adjustments.

## Pediatric Print Arm Sling



Size	Catalog No.
XXSmall	28-4002
XSmall	28-4003
Small	28-4004

### Construction

Pediatric print design.

## Plaid Arm Sling



Size	Catalog No.
XSmall	28-4503
Small	28-4504
Medium	28-4505
Large	28-4506
XLarge	28-4507

### Construction

Heavy plaid material with metal hook closure.

### Features

Closure for ease of application and adjustment.



# KENAD ARM SLINGS

Kenad SG Medical, Inc.

## Universal Arm Elevator



### Description

Economy Version w/ Ice Bags  
Economy Version w/o Ice Bags  
Deluxe Version w/Ice Bags  
Deluxe Version w/o Ice Bags

### Catalog No.

28-5100  
28-5110  
28-5200  
28-5210

### Construction

Beige cotton duck material. Deluxe style is made of heavy blue denim.

### Features

Proven design permits use as either an arm sling or as an arm elevator. Both styles available with or without ice bag.

## Universal Foam Arm Elevator



### Size

Child  
Adult

### Catalog No.

28-5304  
28-5306

### Construction

1" open-cell foam with 1 1/2" cotton webbing.

### Features

Hook and loop closure facilitates application and removal.

## Cradle Arm Sling



### Size

XSmall  
Small  
Medium  
Large

### Catalog No.

28-5403  
28-5404  
28-5405  
28-5406

### Construction

1" wide tie straps. Navy blue polyester/cotton fabric.

### Features

This versatile trapezoidal sling is open at both ends and supports either forearm in a natural, restful position. Tie straps can be applied in three positions: under sling, over sling, or through sling. Washable. No neck pad.

## Universal Hemiplegic Sling



### Description

Universal Hemiplegic Sling

### Catalog No.

28-5500

### Construction

Polyester/cotton material.

### Features

A sturdy, completely adjustable sling designed for the hemiplegic patient. Provides effective support for patients with cerebral vascular accident wounds or for stroke patients. Once sling is properly adjusted to patient for correct forearm angle, it may be applied or removed without assistance. Fits either arm. One size fits all.

# KENAD ARM SLINGS

Kenad SG Medical, Inc.

## Mesh Arm Sling



Size	Catalog No.
Small	28-6004
Medium	28-6005
Large	28-6006
XLarge	28-6007

### Construction

Blue nylon mesh material.

### Features

Hook and loop closure with side buckle to facilitate adjustment. Mesh allows air flow through sling for patient comfort.

## Universal Collar/Cuff



Size	Catalog No.
Adult	28-6500
Child	28-6501

### Construction

1 1/2" cotton/poly webbing and 3/8" foam. Slide buckle for easy adjustment.

### Features

Slide buckle at wrist and neck allows for quick adjustments. Ideal for use in patients who have tendon contractures of the elbow due to long-term cast wear.

## Burbon Sling



Size	Catalog No.
Small	28-7004
Medium	28-7005
Large	28-7006
XLarge	28-7007

### Construction

Heavy blue cotton duck material with 1" webbing. Metal hook closure.

### Features

This design features 1" webbing that passes through two 1" links at mid-sling. Webbing then passes under armpit and over the shoulder. Helps keep shoulder semi-immobilized. Closure allows easy adjustment.



# KENAD ARM SLINGS

Kenad SG Medical, Inc.

## Blue Sta-Level Arm Sling



Size	Catalog No.
XSmall	28-7403
Small	28-7404
Medium	28-7405
Large	28-7406
XLarge	28-7407

### Construction

Polyester/cotton twill. 1 1/2" poly/cotton webbing with hook and loop closure.

### Features

Arm weight is more evenly distributed with wider shoulder strap. Hook and loop closure with slide buckle permits a wide range of adjustments. Sling opens completely for easy insertion of arm.

## Blue Sling w/Neck Pad



Size	Catalog No.
Small	28-7604
Medium	28-7605
Large	28-7606

### Construction

Heavy 10 1/2 oz. royal blue twill. 1 1/2" heavy poly/cotton webbing.

### Features

Foam neck pad for patient comfort. Slide buckle on strap for easy adjustment. Thumb loop allows support of hand. High quality heavy material for durability and patient comfort. Hook and loop closure.

## Blue Arm Sling w/Pilot Hook Closure



Size	Catalog No.
XSmall	28-7503
Small	28-7504
Medium	28-7505
Large	28-7506
XLarge	28-7507

### Construction

1" cotton webbing.

### Features

Metal hook closure provides for easy adjustment. Foam neck pad for patient comfort.

## Khaki Arm Sling



Size	Catalog No.
XSmall	28-7903
Small	28-7904
Medium	28-7905
Large	28-7906
XLarge	28-7907

### Construction

Heavy khaki material.

### Features

A sturdy, completely adjustable sling with hook and loop closure with side buckle.



# KENAD ARM SLINGS

Kenad SG Medical, Inc.

## Denim Arm Sling w/Metal Hook Closure



Size	Catalog No.
XSmall	28-7703
Small	28-7704
Medium	28-7705
Large	28-7706
XLarge	28-7707

### Construction

Heavy denim material.

### Features

Metal hook closure for easy adjustment.

## Adjust/Ease Arm Sling



Size	Catalog No.
XSmall	28-9003
Small	28-9004
Medium	28-9005
Large	28-9006
XLarge	28-9007

### Construction

Heavy 10<sup>1</sup>/<sub>2</sub> oz. navy blue twill. 1<sup>1</sup>/<sub>2</sub>" heavy poly/cotton webbing.

### Features

Easy adjustment with metal hook closure. Foam neck pad for patient comfort. Thumb loop allows support of hand. High quality heavy material for durability.

# KENAD TRACTION PRODUCTS

Kenad SG Medical, Inc.

## Universal Pelvic Traction Belt



**Catalog No.**  
29-1000

### Construction

$\frac{3}{8}$ " foam with hook compatible brushed pile nylon. Hook and loop closure.

### Features

Traction straps may be positioned in various locations on the belt to provide effective, comfortable traction support. The two traction are removable, making the belt either single or double pull style. Universal design reduces inventory requirements.

## Universal Bucks Traction Boot



**Size**  
Pediatric  
Adult

**Catalog No.**  
29-2800  
29-3000

### Construction

Porous eggcrate foam. Aluminum stays with hook and loop closure.

### Features

Designed for use with Russell, Bucks, and other skin applications to the leg. Ten feet of traction cord is included. Universal design reduces inventory requirements.

## Bucks Traction Boot Sized



**Catalog No.**  
29-3103  
29-3104  
29-3105  
29-3106  
29-3107

### Construction

Porous eggcrate foam. Aluminum stays with hook and loop closure.

### Features

Designed for use with Russell, Bucks, and other skin applications to the leg. Ten feet of traction cord is included.

## Complete Cervical Traction Sets



### Description

Cervical Traction Set, Over the Door, Universal  
Cervical Traction Set, Small  
Cervical Traction Set, Medium  
Cervical Traction Set, Large  
Cervical Traction Set, No Halter  
Cervical Traction Set, Dot Head Halter  
Cervical Traction Set, Under Mattress

### Catalog No.

29-3600  
29-3604  
29-3605  
29-3606  
29-3610  
29-3620  
29-3700

### Features

All components necessary for cervical traction from any door or under mattress. Consists of head halter, spreader bar, water weight bag, over-door bracket with pulleys and rope. Head halter size determines the size of the traction set. Universal and Under-Mattress sets include Diskard Head Halter.

# KENAD TRACTION PRODUCTS

Kenad SG Medical, Inc.

## Water Weight Bag



**Description**  
Water  
Weight  
Bag

**Catalog No.**  
29-3630


### Construction

Heavy-duty vinyl.

### Features

Clear, easy-to-read calibrations for determination of correct weight up to 20 pounds.

## Harris Splint Covers



**Description**  
Harris  
Splint  
Covers

**Catalog No.**  
29-3800

### Construction

Plush Kodel.

### Features

Used to cover Harris splints for balanced traction. Sold in sets only.

## Head Halter



**Size**  
Small  
Medium  
Large

**Catalog No.**  
29-3904  
29-3905  
29-3906

### Construction

Strong cotton flannel.

### Features

Slide buckles with catch provide one-time adjustment and easy application and removal.

## Diskard Disposable Head Halter



**Size**  
Universal

**Catalog No.**  
29-4000

### Construction

Surgical gauze padded with polyethylene cellulose.

### Features

Disposable head halter will not irritate chin. Padded to avoid perspiration saturation. Sold in dozen lots only.



# KENAD TRACTION PRODUCTS

Kenad SG Medical, Inc.

## Universal Head Halter



**Description**  
Repro Type

**Catalog No.**  
29-4100

### Construction

$\frac{3}{8}$ " velfoam covered with royal blue cotton duck material outside and soft cotton flannel inside.

### Features

Slide buckle with catch provides one-time adjustment and facilitates application and removal.

## Universal Vinyl Head Halter



**Catalog No.**  
29-4200

### Construction

Brown vinyl laminated to  $\frac{3}{8}$ " foam with tricot lining. Hook and loop closure.

### Features

Designed for use with cervical traction. Universal size reduces inventory requirements.

## Universal Dot Head Halter



**Catalog No.**  
29-4210

### Construction

White canvas with soft cotton backing for comfort. Hook and loop closure.

### Features

Designed for use with cervical traction. Universal size reduces inventory requirements.

## Kodel® Knee Slings



**Size**

Universal	29-4600
Small	29-4704
Medium	29-4705
Large	29-4706
XLarge	29-4707

### Construction

Kodel material with two metal grommets provided for attachment to spreader bar or traction cord.


### Features

Designed for use in Russell traction and other suitable applications.

# KENAD TRACTION PRODUCTS

Kenad SG Medical, Inc.

## Traction Cord



**Length**  
400 yds.  
100 ft.

**Cord  
Diameter**  
 $\frac{5}{32}$ "  
 $\frac{3}{16}$ "

**Material**  
Nylon  
Nylon

**Catalog No.**  
29-5500  
29-5600

### Construction

Available in nylon.

### Features

400-yard nylon roll available in  $\frac{5}{32}$ " diameter only.

## Weil Pelvic Sling



**Description**  
Weil  
Pelvic  
Sling

**Catalog No.**  
29-5800

### Construction

Heavy beige cotton duck material with Kodel lining.

### Features

Four adjustment slots permit proper fitting. used in the treatment of pelvic fractures and for immobilization and stabilization of low back pathology.

# PATIENT PROTECTIVE PRODUCTS

Kenad SG Medical, Inc.

## Protective Mitt



Size	Catalog No.
Adult	50-1000
Child	50-1010
Infant	50-1020
Adult	50-1000BX
Child	50-1010BX
Infant	50-1020BX

### Construction

$\frac{3}{8}$ " white velfoam with hook closure and slide buckle at wrist.

### Features

Prevents patient from grasping intravenous and nasogastric tubes, or from picking and scratching wounds and dressings. Wrist strap with quick-release buckle may be attached to link on mitt to limit movement of limb.

## Straps For Mitts & Vests



### Description

Protective Mitt Wrist Strap	50-1099
Protective Mitt Wrist Strap 10 PR	50-1009BX
Protective Vest Shoulder Strap	50-4099
Protective Vest Shoulder Strap 10 PR	50-4099BX

### Catalog No.

### Construction

1" poly webbing with buckle (male/female).

### Features

For use with all vests and mitts. Vest use is designed to keep patient flat after surgical procedures. Mitt use is designed to limit movement of arms.

## Protective Mitt w/Stays



Size	Catalog No.
Adult	50-1100
Child	50-1110
Infant	50-1120
Adult	50-1100BX
Child	50-1110BX
Infant	50-1120BX

### Construction

$\frac{3}{8}$ " white velfoam with hook closure and slide buckle at wrist. Aluminum stays

### Features

Stays prevent patient from bending hand while inside mitt. Prevents patient from grasping intravenous and nasogastric tubes, or from picking and scratching wounds and dressings. Wrist strap with quick-release buckle may be attached to link on mitt to limit movement of limb.

## White Limb Holder\*



Size	Description	Catalog No.
Universal	White Limb Holder	50-1400
	Threaded White Tails	50-1400BX

### Construction

Soft white foam cuffs. No-slip side buckle wrist closure.

### Features

Buckle wrist closure facilitates cuff application and prevents tourniquet effect. Bed straps attach to bed frame without tying and attach to wrist straps with Delrin quick-release buckle system, safely out of patient's reach


\*U.S. Patent No. 4,422,455



# PATIENT PROTECTIVE PRODUCTS

Kenad SG Medical, Inc.

## Blue Limb Holder\*



Size	Description	Catalog No.
Universal	Blue Limb Holder	50-1500
	Blue Limb Holder Box	50-1500BX

### Construction

Soft blue foam cuffs. No-slip side buckle wrist closure.

### Features

Buckle wrist closure facilitates cuff application and prevents tourniquet effect. Bed straps attach to bed frame without tying and attach to wrist straps with Delrin quick-release buckle system, safely out of patient's reach.

*\*U.S. Patent No. 4,422,455*

## Threaded Loop End Tails



Catalog No.
50-1499
50-1499BX

### Construction

1 1/2" poly webbing with quick release buckle.

### Features

Replacement tails for product number 50-1400 White Limb Holder.

## Loop End White Tails



Description	Catalog No.
Loop End White Tails	50-1599
	50-1599BX

### Construction

1 1/2" poly webbing, female buckle.

### Features

Replacement. For use with product number 50-1500.

## Limb Holder w/Removable Cuff



Size	Description	Catalog No.
Universal	Black Limb Holder	50-1600
Universal	Removable cuff	50-1600BX

### Construction

Foam cuffs attach to wrist strap with Velcro closure.

### Features

Removable cuffs

# PATIENT PROTECTIVE PRODUCTS

Kenad SG Medical, Inc.

## Removable Cuffs


**Description**

Removable Cuffs

**Catalog No.**

50-1655

**Construction**

$\frac{3}{8}$ " foam with hook compatible brushed nylon.  
Hook closure.

**Features**

Replacement cuffs for product. 50-1600.

## Loop End Black Tails


**Description**

Loop End Black Tails  
Loop End Black Tails

**Catalog No.**

50-1699  
50-1699BX

**Construction**

$1\frac{1}{2}$ " poly webbing, female buckle.

**Features**

Replacement. For use with product number  
50-1600.

## Deluxe Pediatric Limb Holder


**Size**

Universal  
Universal

**Description**

Deluxe Pediatric Limb Holder  
Deluxe Pediatric Limb Holder 10 PR

**Catalog No.**

50-1700  
50-1700BX

## Pediatric Loop End Tails


**Catalog No.**

50-1799  
50-1799BX

**Construction**

1" poly webbing with quick release buckle.

**Features**

Replacement tails for use with product number 50-1700 Pediatric Limb Holder.

# PATIENT PROTECTIVE PRODUCTS

Kenad SG Medical, Inc.

## Waist Holder



Size	Description	Catalog No.
Universal	Pediatric Waist Holder	50-3400
Universal	Pediatric Waist Holder 10 EA	50-3400BX
Universal	Adult Waist Holder	50-3500
Universal	Adult Waist Holder 10 EA	50-3500BX
Universal	Adult Waist Holder w/Slit	50-3600
Universal	Adult Waist Holder w/Slit 10 EA	50-3600BX
Universal	Adult Waist Holder X-Wide	50-3700
Universal	Adult Waist Holder X-Wide 10 EA	50-3700BX

### Construction

$\frac{3}{8}$ " foam with hook compatible brushed nylon.  $1\frac{1}{2}$ " poly/cotton webbing and quick-release buckle system.

### Features

Offers comfortable control and safety for bed and wheelchair applications.

## Blue Vest Restraint – No Sleeves



Size	Description	Catalog No.	Catalog No.
XSmall	No Sleeves	50-4003	50-4003BX
Small	No Sleeves	50-4004	50-4004BX
Medium	No Sleeves	50-4005	50-4005BX
Large	No Sleeves	50-4006	50-4006BX
XLarge	No Sleeves	50-4007	50-4007BX
XXLarge	No Sleeves	50-4008	50-4008BX

### Construction

Royal blue 10 oz. cotton twill,  $1\frac{1}{2}$ " poly/cotton webbing and hook and loop closure at waist. Nylon quick-release buckle system.

### Features

In wheelchair application, vest provides posture control for patients who have difficulty sitting erect. Also gives comfortable control with optional shoulder straps in bed application. Quick release buckle system for quick removal and safety. Hook and loop closure at waist helps prevent tourniquet effect.

## Mesh Vest Restraint – No Sleeves



Size	Catalog No.	Catalog No.
XSmall	50-4103	50-4103BX
Small	50-4104	50-4104BX
Medium	50-4105	50-4105BX
Large	50-4106	50-4106BX
XLarge	50-4107	50-4107BX
XXLarge	50-4108	50-4108BX

### Construction

Royal blue 10 oz. mesh,  $1\frac{1}{2}$ " poly/cotton webbing and hook and loop closure at waist. Nylon quick-release buckle system.

### Features

In wheelchair application, vest provides posture control for patients who have difficulty sitting erect. Also gives comfortable control with optional shoulder straps in bed application. Quick-release buckle system for quick removal and safety. Hook and loop closure at waist helps prevent tourniquet effect.



# PATIENT PROTECTIVE PRODUCTS

Kenad SG Medical, Inc.

## Blue Vest Restraint – w/Sleeves



Size	Description	Catalog No.	Catalog No.
XSmall	Sleeved Vest	50-4203	50-4203BX
Small	Sleeved Vest	50-4204	50-4204BX
Medium	Sleeved Vest	50-4205	50-4205BX
Large	Sleeved Vest	50-4206	50-4206BX
XLarge	Sleeved Vest	50-4207	50-4207BX
XXLarge	Sleeved Vest	50-4208	50-4208BX

### Construction

Royal blue 10 oz. twill, 1 1/2" poly/cotton webbing and hook and loop closure at waist. Nylon quick-release buckle system.

### Features

In wheelchair application, vest provides posture control for patients who have difficulty sitting erect. Also gives comfortable control with optional shoulder straps in bed application. Quick-release buckle system for quick removal and safety. Hook and loop closure at waist helps prevent tourniquet effect.

## Mesh Vest Restraint – w/Sleeves



Size	Catalog No.	Catalog No.
XSmall	50-4303	50-4303BX
Small	50-4304	50-4304BX
Medium	50-4305	50-4305BX
Large	50-4306	50-4306BX
XLarge	50-4307	50-4307BX
XXLarge	50-4308	50-4308BX

### Construction

Royal blue 10 oz. mesh, 1 1/2" poly/cotton webbing and hook and loop closure at waist. Nylon quick-release buckle system.

### Features

In wheelchair application, vest provides posture control for patients who have difficulty sitting erect. Also gives comfortable control with optional shoulder straps in bed application. Quick-release buckle system for quick removal and safety. Hook and loop closure at waist helps prevent tourniquet effect.

# PATIENT PROTECTIVE PRODUCTS

Kenad SG Medical, Inc.

## Sleeved Vest Restraint – w/Tails



Size	Description	Catalog No.	Catalog No.
XSmall	Sleeved Vest w/Tails	50-4403	50-4403BX
Small	Sleeved Vest w/Tails	50-4404	50-4404BX
Medium	Sleeved Vest w/Tails	50-4405	50-4405BX
Large	Sleeved Vest w/Tails	50-4406	50-4406BX
XLarge	Sleeved Vest w/Tails	50-4407	50-4407BX
XXLarge	Sleeved Vest w/Tails	50-4408	50-4408BX

### Construction

Royal blue 10 oz. twill, 1½" poly/cotton webbing and hook and loop closure at waist. Nylon quick-release buckle system.

### Features

In wheelchair application, vest provides posture control for patients who have difficulty sitting erect. Also gives comfortable control with optional shoulder straps in bed application. Quick-release buckle system for quick removal and safety. Hook and loop closure at waist helps prevent tourniquet effect.

## Pediatric Vest Restraint



Size	Description	Catalog No.	Catalog No.
Small	Pediatric Vest	50-4604	50-4604BX
Medium	Pediatric Vest	50-4605	50-4605BX
Large	Pediatric Vest	50-4606	50-4606BX

### Construction

Royal blue 10 oz. twill, 1½" poly/cotton webbing and hook and loop closure at waist. Nylon quick-release buckle system.

### Features

In wheelchair application, vest provides posture control for patients who have difficulty sitting erect. Also gives comfortable control with optional shoulder straps in bed application. Quick-release buckle system for quick removal and safety. Hook and loop closure at waist helps prevent tourniquet effect.

## Elbow & Wrist Immobilizer



Size	Catalog No.	Catalog No.	Dimension
XXSmall	50-5002	50-5002BX	4" x 7"
XSmall	50-5003	50-5003BX	5" x 5"
Small	50-5004	50-5004BX	7" x 13"
Medium	50-5005	50-5005BX	8" x 15"
Large	50-5006	50-5006BX	10" x 18"
XLarge	50-5007	50-5007BX	13" x 18"

### Construction

¾" white velfoam with aluminum stays positioned vertically.

### Features

One of our most popular products. Stays prevent bending of the wrist and elbow to keep patient from pulling out upper extremity tubes. Can also be used as an I.V. board by cutting a hole in the foam and wrapping around I.V. tube. Hole allows monitoring for infiltration of I.V.



# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Telemetry Holder



**Size**  
Universal

**Catalog No.**  
61-1010

### Construction

$\frac{3}{8}$ " white velfoam with hook tab closure.

### Features

Great for holding telemetry units to prevent unit from falling from bed or disconnecting by patient. Designed to hold telemetry holder and prevent unit from damage.

## Tracheotomy Tube Holder



**Description**  
Tracheotomy Tube Holder

**Catalog No.**  
61-1100

### Construction

$\frac{1}{8}$ " velfoam with 12" strip adhesive mole skin inside.

### Features

Excellent for holding endo tubes in place.

## Kodel® Elbow Pad



**Description**  
Kodel® Elbow Pad

**Catalog No.**  
61-2100

### Features

A soft, washable plush Kodel® elbow protector contoured for patient comfort. Designed to facilitate air circulation and to distribute even pressure on the elbow.

## Convuluted Foam Elbow & Heel Pad



**Description**  
Convuluted Foam Elbow & Heel Pad

**Catalog No.**  
61-2200

### Construction

Medical grade polyurethane foam.

### Features

Designed to provide comfort and to discourage decubitus ulcers for long-term immobile patients. Convoluted surface allows greater air circulation.



# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Tube Tac Tube Positioning System



### Description

Universal Tube Tac System

### Catalog No.

61-2500

### Features

Positions and secures catheters, nasal, and IV tubes to limbs or head. Holds tube firmly with adhesive surface. Hook fasteners grasps pile of positioning strap.

## Gait Belt w/Buckle Closure



### Size

Universal 4" x 4"  
Universal 4" x 5"

### Catalog No.

61-2710  
61-2712

### Construction

$\frac{3}{8}$ " foam laminated to soft hook compatible brushed nylon. 1 $\frac{1}{2}$ " quick-release buckle. Hook and loop closure.

### Features

Designed to be used to assist in ambulation or as in aid in patient transfer.

## Abduction Pillow – Leg



### Style

Infant  
Youth  
Child  
Adult Small  
Adult Medium  
Adult Large

### Catalog No.

61-3001  
61-3002  
61-3003  
61-3004  
61-3005  
61-3006


### Construction

Sturdy polyurethane foam.

### Features

Maintains hip in abduction. Supports and aligns patient's leg. Foam may be trimmed to modify for bony prominences.

## Universal Hoyer Sling



### Size

Universal  
Hoyer Sling

### Catalog No.

61-3200

### Construction

Blue nylon mesh reinforced with 1 $\frac{1}{2}$ " cotton webbing.

### Features

Replacement sling for Hoyer lift.

# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Foot Block Elevator



**Dimensions**  
8" x 8" x 12"

**Catalog No.**  
61-4200

### Construction

Open cell polyurethane foam with hook and pile strap.

### Features

Slit allows easy insertion of limb.

## Kodel® Heel Protectors



**Style**  
Deluxe  
Economy

**Catalog No.**  
61-4600  
61-4610

### Construction

Plush Kodel® with hook and loop closure.

### Features

Available in deluxe and economy versions. Designed to discourage decubitus ulcers for long-term immobile patients. Universal design reduces inventory requirements. Sold in pairs only.

## Heel & Elbow Pads



**Size**  
Small  
Medium  
Large  
XLarge

**Catalog No.**  
61-4704  
61-4705  
61-4706  
61-4707

### Construction

Tubular elastic, foam pad.

### Features

Provides comfort and discourages decubitus ulcers for long-term immobile patients. Distributes even pressure on the heel or elbow.

## Heavy Duty Foot Cradle



**Description**  
Heavy Duty  
Foot Cradle

**Catalog No.**  
61-4800

### Construction

Convuluted foam laminated to 3/8" blue velfoam.


### Features

Designed for long-term decubitus prevention for the feet. Sold in singles only.

# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Bed Rail Pads



### Description

Bed Rail Pads  
Sides Only  
Head (1)

### Size

Child  
Adult  
Adult

### Catalog No.

61-5500  
61-5600  
61-5610


### Construction

1" medium density foam covered with cotton stockinette.

### Features

Hook and loop straps allow secure fit to side rails of head board.

## Replacement CPM Pads – Sutter



### Description

Sutter Pads

### Catalog No.

61-6200

### Construction

Plush Kodel® with hook and loop closures.

### Features

Replacement pads for Sutter Continuous Passive Motion machines. Sold in sets only.

## Replacement CPM Pads – Richards



### Description

Richards  
Pads

### Catalog No.

61-6210

### Construction

Plush Kodel® with hook and loop closures.

### Features

Replacement pads for Richards Continuous Passive Motion machines. Sold in sets only.

## Replacement CPM Pads – Danninger Model 300



### Description

Danninger  
Model 300  
Pads

### Catalog No.

61-6220

### Construction

Plush Kodel® with hook and loop closures.

### Features

Replacement pads for Danninger Model 300 Continuous Passive Motion machines. Sold in sets only.



# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Stryker CPM Pads



### Description

Universal  
7 Pieces  
4 Pieces

### Catalog No.

61-6240  
61-6260

### Construction

Plush Kodel® with hook and loop closures.

### Features

Replacement pads for Stryker Continuous Passive Motion machines. Sold in sets only.

## CPM Pads Stryker / Danninger



### Description

CPM Pads

### Catalog No.

61-6250

### Construction

Plush Kodel. Hook and loop closure.

### Features

Replacement pad for continuous passive motion machines. Designed to fit Danninger or Stryker model.

## Powerflex CPM Pads



### Description

Powerflex  
CPM Pads

### Catalog No.

61-6270

### Construction

Plush Kodel® with hook and loop closures.

### Features

Continuous Passive Motion machine replacement pads. Sold in sets only.

## Danninger Model 400 CPM Pads



### Description

Danninger  
Model 400  
Pads

### Catalog No.

61-6280

### Construction

Plush Kodel® with hook and loop closures.

### Features

Replacement pads for Danninger Model 300 Continuous Passive Motion machines. Sold in sets only.

# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Optiflex CPM Pad



**Description**  
CPM Pads

**Catalog No.**  
61-6285

### Construction

Plush Kodel. Hook and loop closure and  $\frac{3}{8}$ " foam laminate.

### Features

Replacement pad for continuous passive motion machines. Designed to fit Optiflex model.

## Thomas Splint / CPM Pads



**Description**  
CPM Pads

**Catalog No.**  
61-6290

### Construction

Plush Kodel. Hook and loop closure and  $\frac{3}{8}$ " foam laminate.

### Features

Replacement pad for continuous passive motion machines. Designed to fit Thomas model.

## Kodel® Bed Pad



### Size

36" x 40" 61-6600  
24" x 30" 61-6601  
30" x 40" 61-6602  
30" x 60" 61-6603  
36" x 80" 61-6700

### Catalog No.

### Construction

Soft Kodel®.

### Features

Bed pad helps reduce the friction and pressure that contribute to the formation of decubitus ulcers.

## Foam Block Arm Elevator



**Description**  
Foam Block  
Arm Elevator

**Catalog No.**  
61-6900

### Construction

Open cell polyurethane foam.

### Features

Made with slit on one side for insertion of arm. Holes on side permit air circulation.

# PATIENT ASSIST PRODUCTS

Kenad SG Medical, Inc.

## Kenad Spinal Frame



Description	Dimensions	Size	Catalog No.
Frame	Universal	Universal	61-9900
Bolster	10 <sup>1</sup> / <sub>2</sub> " x 20"	Large	61-9910
Bolster	10 <sup>1</sup> / <sub>2</sub> " x 18"	Medium	61-9920
Bolster	10 <sup>1</sup> / <sub>2</sub> " x 16"	Small	61-9930
Bolster	7 <sup>1</sup> / <sub>2</sub> " x 20"	Large	61-9940
Bolster	7 <sup>1</sup> / <sub>2</sub> " x 18"	Medium	61-9950
Bolster	7 <sup>1</sup> / <sub>2</sub> " x 16"	Small	61-9960

### Construction

Wooden frame. Foam.

### Features


Frame accommodates x-ray cassette. Adjustable foam bolsters aid in patient positioning.



# SPLINTS

Kenad SG Medical, Inc.

## Nose Splint



**Description**  
Nose  
Splint

**Catalog No.**  
62-0800

### Construction

Padded aluminum.

### Features

Nose splint is malleable, can be formed to contour.

## Finger/Thumb Splint



**Dimensions**  
1" x 7"  
1" x 9"  
1" x 12"

**Size**  
Small  
Medium  
Large

**Catalog No.**  
62-0904  
62-0905  
62-0906

### Construction

Vinyl coated aluminum.

### Features

May be shaped to maintain finger or thumb in extreme flexion or extension. For maximum adjustability, the center can be moved up and down the entire length of the splint.

## Toad Splints



**Length**  
2 1/4"  
2 3/4"  
2 7/8"

**Size**  
Small  
Medium  
Large

**Catalog No.**  
62-1004  
62-1005  
62-1006

### Construction

Malleable aluminum. 1/4" thick polyester padding.

### Features

Aluminum arms fold to hold finger in desired alignment. Maintains interphalangeal position for injuries to distal joints. Padding keeps finger dry and comfortable. Padded style only. Packaged in quantities of 12.

## Padded Finger Splint Kit



### Kit includes

12 – Cat. No. 62-1006 – Large Toad Splints  
12 – Cat. No. 62-1005 – Medium Toad Splints  
12 – Cat. No. 62-1302 – Padded 6" Finger Splints  
6 – Cat. No. 62-1303 – Padded 2 1/2" Finger Splints  
6 – Cat. No. 62-1304 – Padded 1 1/2" Finger Splints  
9 – Cat. No. 62-1704 – Padded 3 1/2" Finger Cots  
6 – Cat. No. 62-1703 – Padded 3" Finger Cots  
6 – Cat. No. 62-1702 – Padded 2" Finger Cots

### Catalog No.

62-1300

### Features

An assortment of 69 padded finger splints packaged in an attractive storage container. When closed, this kit becomes a handy carrying case. Ideal for doctor's office, clinics, emergency rooms, and industrial first aid departments.

# SPLINTS

Kenad SG Medical, Inc.

## Padded Finger Splint



	7/8" Width Padded Catalog No.	7/8" Width Unpadded Catalog No.	5/8" Width Padded Catalog No.	5/8" Width Unpadded Catalog No.
Length				
7"	62-1301P	62-1301	62-1501P	62-1501
6"	62-1302P	62-1302	62-1502P	62-1502
2 1/2"	62-1303P	62-1303	62-1503P	62-1503
1 1/2"	62-1304P	62-1304	62-1504P	62-1504
4"	62-1305P	62-1305	62-1505P	62-1505

### Construction

1/4" polyester foam bonded to aluminum.

### Features

Curved to conform to the finger in a semi-flexed position. Used in transverse fractures of the phalanges where extension is not required. Packaged in quantities of 12.

## Finger Cot



Size	Padded Catalog No.	Unpadded Catalog No.
1"	62-1701P	62-1701
2"	62-1702P	62-1702
3"	62-1703P	62-1703
3 1/2"	62-1704P	62-1704

### Construction

1/4" thick polyester foam bonded to aluminum.

### Features

Used for protecting injured finger tips, especially after amputation, thus eliminating heavy dressings and affording more adequate ventilation to injured finger. Packaged in quantities of 12.

## Finger Protectors



Size	Padded Catalog No.	Unpadded Catalog No.
Small	62-1904P	62-1904
Medium	62-1905P	62-1905
Large	62-1906P	62-1906

### Features

Provides end protection to injured fingers and eliminates bulky dressing. Packaged in quantities of 12.



# SPLINTS

Kenod SG Medical, Inc.

## Foam Padded Finger Strips



Size	Catalog No.	Qty./Pkg.
1" x 18"	62-2101	12
3/4" x 18"	62-2102	12
1/2" x 18"	62-2103	12
1/2" x 9"	62-2104	24
3/4" x 9"	62-2105	24
1" x 9"	62-2106	24
1 1/2" x 18"	62-2107	12
2" x 18"	62-2108	12
3" x 18"	62-2109	12

### Features

Used for fractures of the phalanges and may be applied as a baseball splint, finger protector, or finger cot. Easy to form and can be cut to desired length with scissors. Sold in full packages only.

## Stax Finger Splints



Size	Catalog No.	Qty./Pkg.
1	62-2501	12
2	62-2502	12
3	62-2503	12
4	62-2504	12
5	62-2505	12
5 1/2"	62-2506	12
6	62-2507	12
7	62-2508	12
Strap	62-2509	25/Box

### Construction

Skin colored plastic.

### Features

Designed to support the distal joint of finger in extension permitting unrestricted movement of the proximal interphalangeal joint. The Stax splint is normally used without padding and is fixed to the middle phalanx with adhesive strip, automatically elevating the distal joint to its proper extended position. In cases where the finger is too tender, a large size splint with padding on the inside should be used until it can be replaced with a normal size Stax splint.

## Stax Finger Splint Kit



### Kit includes

4 ea. - Cat. No. 62-2501  
4 ea. - Cat. No. 62-2502  
4 ea. - Cat. No. 62-2503  
5 ea. - Cat. No. 62-2504  
4 ea. - Cat. No. 62-2505  
4 ea. - Cat. No. 62-2506  
4 ea. - Cat. No. 62-2507  
1 ea. - Cat. No. 62-1702

**Catalog No.**  
62-2300

### Features

An assortment of 30 Stax splints packaged in an attractive storage container.

## Malleable Baseball Finger Splint



Size	Catalog No.
Small	62-2704
Medium	62-2705
Large	62-2706

### Construction

Padded styles. Aluminum.

### Features

May be used for fractures of the distal phalanx. Easily molded to maintain desired flexion of fractured phalanx.



# SPLINTS

Kenad SG Medical, Inc.

## Cockup Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult	63-1000P	63-1000
Youth	63-1010P	63-1010
Child	63-1020P	63-1020

### Construction

Vinyl coated aluminum.

### Features

The ball shape in the hand affords a comfortable rest and maintains the arch of the palm. The degree of dorsiflexion may be changed by bending at the indicated point.

## Colles Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult Right	63-1201P	63-1201
Adult Left	63-1202P	63-1202
Youth Right	63-1203P	63-1203
Youth Left	63-1204P	63-1204
Child Right	63-1205P	63-1205
Child Left	63-1206P	63-1206

### Construction

Vinyl coated aluminum.

### Features

Generally used for distal end fractures of the ulna and radius. Can be fitted to correct angle of flexion by slightly twisting and bending at the indicated locations. Hand-piece is solid aluminum contoured to fit the hand.

## Special Colles Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult Right	63-1401P	63-1401
Adult Left	63-1402P	63-1402
Youth Right	63-1403P	63-1403
Youth Left	63-1404P	63-1404
Child Right	63-1405P	63-1405
Child Left	63-1406P	63-1406
Infant Right	63-1407P	63-1407
Infant Left	63-1408P	63-1408

### Construction

Vinyl coated aluminum.

### Features

This splint conforms to the shape of the arm and hand in a supinated and semipronated posture. The hand is maintained in the palmar flexed and abducted position.

## Forearm Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult Right	63-1601P	63-1601
Adult Left	63-1602P	63-1602
Youth Right	63-1603P	63-1603
Youth Left	63-1604P	63-1604
Child Right	63-1605P	63-1605
Child Left	63-1606P	63-1606
Infant Right	63-1607P	63-1607
Infant Left	63-1608P	63-1608

### Construction

Vinyl coated aluminum.

### Features

Useful in fractures of the ulna and radius. It can also be used for synovitis of wrist. Also for sprains and for Colles fractures where abducted position of hand is not desired. It can be molded and shaped to forearm contour.

# SPLINTS

Kenod SG Medical, Inc.

## Metacarpal Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult Right	63-1801P	63-1801
Adult Left	63-1802P	63-1802
Youth Right	63-1803P	63-1803
Youth Left	63-1804P	63-1804
Child Right	63-1805P	63-1805
Child Left	63-1806P	63-1806
Infant Right	63-1807P	63-1807
Infant Left	63-1808P	63-1808

### Construction

Vinyl coated aluminum.

### Features

This splint is useful in the reduction of radial, ulnar, metacarpal, and phalangeal fractures. It assures comfortable immobilization of the hand and forearm.

## Metacarpal Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult Right	63-2001P	63-2001
Adult Left	63-2002P	63-2002
Youth Right	63-2003P	63-2003
Youth Left	63-2004P	63-2004
Child Right	63-2005P	63-2005
Child Left	63-2006P	63-2006

### Construction

Vinyl coated aluminum.

### Features

Designed to fit the hand. A great amount of traction is possible by affixing tape to and splitting over finger. end of tape is then attached to base of splint. As splint is pressured back into palm, traction is exerted on finger. Additional tape holds splint in place.

## Posterior Elbow Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult	63-2400P	63-2400
Youth	63-2410P	63-2410
Child	63-2420P	63-2420

### Construction

Vinyl coated aluminum.

### Features

Eliminates pressure on the olecranon process when used to immobilize the elbow. Bar may be bent to change the angle. Fits either arm.

## Posterior Knee Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult	63-2600P	63-2600
Youth	63-2610P	63-2610
Child	63-2620P	63-2620

### Construction

Vinyl coated aluminum.

### Features

For treatment of fractures of the patella or whenever immobilization of knee is desired. Provides length and rigidity necessary to hold knee during dressing applications.



# SPLINTS

Kenad SG Medical, Inc.

## Posterior Leg Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult	63-2800P	63-2800
Youth	63-2810P	63-2810
Child	63-2820P	63-2820

### Construction

Vinyl coated aluminum.

### Features

Used for fractures of the femur and for immobilization of the knee. Rigid aluminum construction, yet splint can be formed to the contour of the leg. Hep length.

## Posterior Leg Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult	63-3000P	63-3000
Youth	63-3010P	63-3010
Child	63-3020P	63-3020

### Construction

Vinyl coated aluminum.

### Features

Used in treatment of transverse fractures of the tibia, or in fibula and ankle fractures, especially where extension is required. Also used to immobilize the lower leg during dressing application or while transporting patients with injuries. Knee length.

## Mason-Allen Hand Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Large	63-3200P	63-3200
Small	63-3210P	63-3210

### Construction

Vinyl coated aluminum.

### Features

Provides immobilization of entire hand and wrist after hand surgery, or nerve or tendon injuries.

## Speed Hand Splint



Size	Padded Catalog No.	Unpadded Catalog No.
Adult Right	63-3401P	63-3401
Adult Left	63-3402P	63-3402
Youth Right	63-3403P	63-3403
Youth Left	63-3404P	63-3404
Child Right	63-3405P	63-3405
Child Left	63-3406P	63-3406

### Construction

Vinyl coated aluminum.

### Features

Used in the treatment of injuries where it is important to keep the forearm and thumb immobilized in a natural position.



# REHABILITATION

Kenod SG Medical, Inc.

## Isoton Exercise Office Unit



<b>Description</b>	<b>Catalog No.</b>
Complete Office Unit	70-1000

### Features

The Isoton Exercise System is designed for the prevention and rehabilitation of acute and post-op injuries and repetitive motion disorders. The only system on the market that exercises all aspects of the upper and lower extremities, including the shoulder, elbow, wrist, hand, hip, knee, and ankle. The office unit allows the therapist to instruct the patient for home use.

## Isoton Exercise Home Unit



<b>Description</b>	<b>Catalog No.</b>
Complete Home Unit	70-1010

### Features

The home unit offers a complete upper and lower extremity exercise program for the home for faster recovery time. The home unit comes with handy door bracket, six exercise bands, and the leg and hand strap. Eliminates cumbersome knots and needless wear on door knobs. The home unit allows the patient to work at his or her own pace, as directed by their doctor or therapist.



## Isoton Accessories



<b>Description</b>	<b>Catalog No.</b>
Universal Leg & Hand Strap	70-1020
Exercise Band Yellow	70-1030
Exercise Band Red	70-1040
Exercise Band Green	70-1050
Exercise Band Blue	70-1060
Exercise Band Black	70-1070
Exercise Band Silver	70-1080
Door Door Bracket	70-1090
Isoton Exercise Brochure	70-2000
Isoton Exercise Manual	70-2010

### Features

Replacement Parts for the Isoton Office or home unit.

- ▶ Velcro® is a registered trademark of Velcro USA, Inc.
- ▶ The Philadelphia Cervical Collar® is a registered trademark of Charles Grenier and Company, Inc.
- ▶ Kodel® is a registered trademark of Eastman Kodak Company
- ▶ All other Trademarks are the rights of the respected owners